

 HENRY SCHEIN® 
SOLUTIONS FOR HEALTH CARE PROFESSIONALS

Digital DENTISTRY

Your source for **CAD/CAM** professionals



Rely on Us™



 1.800.688.5558

 www.henryschein.ca


TABLE OF CONTENTS

CAD/CAM Accessories.....	45-47
Characterization.....	26-28
Etching / Bonding.....	41
Finishing & Polishing.....	32-35
Furnance Accessories	31
Furnances	30
Implant Support Products	29
Impression Material	45
Milling Blocks	4-25, 48
Permanent Cements	37-40

STAY CONNECTED *with*



Stay informed on the most up-to-date news in the dental industry with our many social network platforms. Just one more reason to **rely on Henry Schein.**



© 2021 Henry Schein Canada. No copying without permission. Not responsible for typographical errors.



Invest in Your Practice...

NOW!

with the Route 66 Lease Program

- No payments for 6 months
- \$99 for the next 6 months
- Then level payments for 48 months*

Henry Schein Financial Services offers:

- Simple application process
- Flexible programs
- Affordable payments
- Fast approval process



To learn more, contact Henry Schein Financial Services



1-800-268-6497 • hsfs@henryschein.ca

*Minimum terms and conditions. Minimum transaction size: \$5,000. All programs feature a \$10 buy out. All transactions are subject to the satisfaction of underwriting guidelines, credit approval by third party lenders, and documentation requirements, and not all applicants will qualify. Certain other restrictions and additional terms and conditions apply. Offer expires December 31, 2021.

©2021 Henry Schein. Neither Henry Schein, Inc. nor Henry Schein Financial Services is a bank and neither conducts banking activities. Henry Schein Financial Services may receive a marketing fee from the vendor for products/services purchased.

BE YOUR BEST

with a NEW end-to-end CAD/CAM workflow
for great restorative treatments



A trusted connection between 3Shape's award-winning TRIOS® innovation and Ivoclar Vivadent's PrograMill® One® now allows you to craft high-end IPS e.max® restorations and give your patients the smiles they deserve.

SCAN
with 3Shape TRIOS® 4



DESIGN
with 3Shape TRIOS
Design Studio®



PRODUCE
great IPS e.max® restorations with
Ivoclar Vivadent PrograMill® One®

Contact your Henry Schein Sales Consultant for more information.
LEARN MORE: 3shape.ivoclarvivadent.com/be-your-best

3shape 


**ivoclar
vivadent**®
passion vision innovation

©2020 Ivoclar Vivadent, Inc. Ivoclar Vivadent, PrograMill, One and IPS e.max are trademarks of Ivoclar Vivadent, Inc. TRIOS and TRIOS Design Studio are trademarks of 3Shape A/S.

IPS e.max® CAD Blocks

Ivoclar Vivadent

High-strength, lithium-disilicate, glass-ceramic blocks. Biaxial flexural strength of 530 MPa. Combine highly aesthetic qualities with exceptional user-friendliness. Provides a wide range of translucency levels, shades, and block sizes to offer greater flexibility for challenging clinical situations. Offers a comprehensive spectrum of indications, which includes thin veneers, veneers, inlays, onlays, partial crowns, crowns, anterior bridges, multiple-unit bridges, 3-unit bridges (up to the second premolar as the terminal abutment), hybrid abutments, and abutment crowns. Furnace required postmilling.



HT (High Translucency)

Ideal for fabricating smaller restorations, inlays, onlays and veneers. Lifelike fluorescence and surround dentition refine the shade of restoration to achieve a balanced true-to-nature chameleon effect.

LT (Low Translucency)

Characterized by their high level of brightness and true-to-nature chroma. Suitable for the fabrication of larger restorations (partial crowns, anterior and posterior crowns and veneers).

MT (Medium Translucency)

Used in cases where a brighter and more translucent look is needed than that imparted by the LT blocks. Restorations made of the MT material are ideal for the staining and cut-back techniques.

MO (Medium Opacity)

For the fabrication of substructures that are placed on vital or slightly discolored prepared teeth. Form an excellent base for lifelike restorations that are completed with the layering technique.

Impulse (Opal)

Provide a combination of life-like opalescence and strength. Ideal for fabricating thin veneers and veneers for light teeth, which require an opalescent effect.

Abutment

With a prefabricated interface. Suitable for the creation of hybrid abutments.

For use with **CEREC®** and **inLab®**



HT Blocks

Size	I12	C14	B40	B40L
Shade				
A1	(945-1112)	(945-1223)	(945-1491)	(945-1504)
A2	(945-1114)	(945-1225)	(945-1495)	(945-1505)
A3	(945-1115)	(945-1226)	(945-1496)	(945-1506)
A3.5	(945-1116)	(945-1230)	(945-1497)	(945-1508)
A4	(945-1117)	(945-1231)	—	—
B1	(945-1118)	(945-1232)	(945-1498)	(945-1510)
B2	(945-1119)	(945-1233)	(945-1499)	(945-1513)
B3	(945-1120)	(945-1234)	—	—
B4	(945-1121)	(945-1235)	—	—
C1	(945-1122)	(945-1236)	(945-1500)	(945-1514)
C2	(945-1123)	(945-1237)	(945-1501)	(945-1515)
C3	(945-1124)	(945-1238)	—	—
C4	(945-1126)	(945-1239)	—	—
D2	(945-1127)	(945-1240)	(945-1502)	(945-1516)
D3	(945-1128)	(945-1242)	—	—
D4	(945-1130)	(945-1243)	—	—
BL1	(945-1131)	(945-1245)	—	—
BL2	(945-1132)	(945-1246)	(945-1503)	(945-1517)
BL3	(945-1133)	(945-1247)	—	—
BL4	(945-1134)	(945-1248)	—	—
	5 per Box		3 per Box	3 per Box

LT Blocks

Size	I12	C14	C16	B32
Shade				
A1	(945-1091)	(945-1198)	(945-1264)	(945-1286)
A2	(945-1092)	(945-1199)	(945-1267)	(945-1287)
A3	(945-1093)	(945-1200)	(945-1268)	(945-1290)
A3.5	(945-1094)	(945-1201)	(945-1270)	(945-1291)
A4	(945-1095)	(945-1203)	—	—
B1	(945-1096)	(945-1204)	(945-1271)	(945-1292)
B2	(945-1097)	(945-1206)	(945-1272)	(945-1293)
B3	(945-1098)	(945-1207)	—	—
B4	(945-1099)	(945-1208)	—	—
C1	(945-1100)	(945-1209)	(945-1275)	(945-1294)
C2	(945-1102)	(945-1210)	(945-1278)	(945-1295)
C3	(945-1103)	(945-1211)	—	—
C4	(945-1104)	(945-1212)	—	—
D2	(945-1105)	(945-1213)	(945-1279)	(945-1300)
D3	(945-1106)	(945-1214)	—	—
D4	(945-1107)	(945-1218)	—	—
BL1	(945-1108)	(945-1219)	—	—
BL2	(945-1109)	(945-1220)	(945-1280)	(945-1301)
BL3	(945-1110)	(945-1221)	—	—
BL4	(945-1111)	(945-1222)	—	—
	5 per Box		5 per Box	3 per Box

IPS e.max® CAD Blocks (Cont'd)

For use with **CEREC®** and **inLab®**

LT Abutment Blocks

Size	A14 L	A14 S	A16 L	A16 S
Shade				
A1	(945-1532)	(945-1518)	(945-1567)	(945-1555)
A2	(945-1533)	(945-1519)	(945-1568)	(945-1556)
A3	(945-1534)	(945-1520)	(945-1569)	(945-1557)
A3.5	(945-1535)	(945-1521)	(945-1570)	(945-1558)
B1	(945-1536)	(945-1523)	(945-1571)	(945-1559)
B2	(945-1538)	(945-1525)	(945-1573)	(945-1560)
C1	(945-1539)	(945-1526)	(945-1574)	(945-1561)
C2	(945-1541)	(945-1528)	(945-1575)	(945-1562)
D2	(945-1542)	(945-1529)	(945-1576)	(945-1563)
BL2	(945-1543)	(945-1531)	(945-1565)	(945-1564)
	5 per Box		5 per Box	

MO Blocks

Size	C14
Shade	
MO0	(945-1634)
MO1	(945-1635)
MO2	(945-1636)
MO3	(945-1637)
MO4	(945-1639)
	5 per Box

MO Abutment Blocks

Size	A14 L	A14 S
Shade		
MO0	(945-1549)	(945-1544)
MO1	(945-1550)	(945-1545)
MO2	(945-1551)	(945-1546)
MO3	(945-1552)	(945-1547)
MO4	(945-1554)	(945-1548)
	5 per Box	

MT Blocks

Size	C14
Shade	
A1	(945-1252)
A2	(945-1253)
A3	(945-1255)
B1	(945-1256)
BL2	(945-1258)
BL3	(945-1260)
BL4	(945-1261)
	5 per Box

Impulse Blocks

Size	C14
Shade	
Opal-O1	(945-1249)
Opal-O2	(945-1251)
	5 per Box

For use with **Planmeca PlanMill®**

HT Blocks

Size	I12	C14
Shade		
A1	(945-0793)	(945-0814)
A2	(945-0794)	(945-0815)
A3	(945-0795)	(945-0817)
A3.5	(945-0796)	(945-0818)
A4	(945-0797)	(945-0819)
B1	(945-0802)	(945-0821)
B2	(945-0803)	(945-0822)
B3	(945-0804)	(945-0823)
B4	(945-0805)	(945-0824)
C1	(945-0806)	(945-0825)
C2	(945-0807)	(945-0826)
C3	(945-0808)	(945-0827)
C4	(945-0809)	(945-0828)
D2	(945-0810)	(945-0829)
D3	(945-0812)	(945-0831)
D4	(945-0813)	(945-0832)
BL1	(945-0833)	(945-0838)
BL2	(945-0834)	(945-0839)
BL3	(945-0836)	(945-0840)
BL4	(945-0837)	(945-0841)
	5 per Box	

LT Blocks

Size	I12	C14	C16	B32
Shade				
A1	(945-0746)	(945-0768)	(945-0859)	(945-0847)
A2	(945-0747)	(945-0770)	(945-0860)	(945-0848)
A3	(945-0748)	(945-0771)	(945-0861)	(945-0849)
A3.5	(945-0749)	(945-0773)	(945-0862)	(945-0850)
A4	(945-0750)	(945-0774)	-	-
B1	(945-0751)	(945-0775)	(945-0864)	(945-0851)
B2	(945-0752)	(945-0776)	(945-0865)	(945-0852)
B3	(945-0753)	(945-0777)	-	-
B4	(945-0754)	(945-0779)	-	-
C1	(945-0755)	(945-0780)	(945-0866)	(945-0853)
C2	(945-0756)	(945-0781)	(945-0869)	(945-0855)
C3	(945-0757)	(945-0782)	-	-
C4	(945-0759)	(945-0783)	-	-
D2	(945-0760)	(945-0784)	(945-0870)	(945-0856)
D3	(945-0761)	(945-0786)	-	-
D4	(945-0762)	(945-0787)	-	-
BL1	(945-0764)	(945-0789)	-	-
BL2	(945-0765)	(945-0790)	(945-0857)	(945-0844)
BL3	(945-0766)	(945-0791)	-	-
BL4	(945-0767)	(945-0792)	-	-
	5 per Box		5 per Box	3 per Box

MT Blocks

Size	C14
Shade	
A1	(945-0643)
A2	(945-0644)
A3	(945-0645)
B1	(945-0646)
BL2	(945-0648)
BL3	(945-0649)
BL4	(945-0650)
	5 per Box

Impulse Blocks

Size	C14
Shade	
Opal-O1	(945-0842)
Opal-O2	(945-0843)
	5 per Box

IPS e.max® CAD Blocks (Cont'd)



For use with **PrograMill®**



HT Blocks

Size	I12	C14
Shade		
A1	(945-1840)	(945-1818)
A2	(945-1841)	(945-1819)
A3	(945-1842)	(945-1820)
A3.5	(945-1843)	(945-1821)
A4	(945-1844)	(945-1822)
B1	(945-1845)	(945-1824)
B2	(945-1846)	(945-1825)
B3	(945-1847)	(945-1826)
B4	(945-1848)	(945-1827)
C1	(945-1851)	(945-1828)
C2	(945-1852)	(945-1830)
C3	(945-1853)	(945-1831)
C4	(945-1854)	(945-1832)
D2	(945-1856)	(945-1833)
D3	(945-1857)	(945-1834)
D4	(945-1864)	(945-1835)
BL1	(945-1836)	(945-1814)
BL2	(945-1837)	(945-1815)
BL3	(945-1838)	(945-1816)
BL4	(945-1839)	(945-1817)
	5 per Box	

LT Blocks

Size	C14	C16
Shade		
A1	(945-1741)	(945-1877)
A2	(945-1742)	(945-1878)
A3	(945-1743)	(945-1879)
A3.5	(945-1744)	(945-1880)
A4	(945-1746)	-
B1	(945-1752)	(945-1881)
B2	(945-1755)	(945-1882)
B3	(945-1760)	-
B4	(945-1770)	-
C1	(945-1771)	(945-1883)
C2	(945-1791)	(945-1884)
C3	(945-1792)	-
C4	(945-1793)	-
D2	(945-1794)	(945-1885)
D3	(945-1803)	-
D4	(945-1805)	-
BL1	(945-1736)	-
BL2	(945-1737)	(945-1876)
BL3	(945-1739)	-
BL4	(945-1740)	-
	5 per Box	

MO Blocks

Size	C14
Shade	
0	(945-1722)
1	(945-1728)
2	(945-1730)
3	(945-1731)
4	(945-1735)
	5 per Box

MT Blocks

Size	C14
Shade	
A1	(945-1809)
A2	(945-1811)
A3	(945-1812)
B1	(945-1813)
BL2	(945-1806)
BL3	(945-1808)
BL4	(945-1807)
	5 per Box

Impulse Blocks

Size	C14
Shade	
Opal-O1	(945-1720)
Opal-O2	(945-1721)
	5 per Box



IPS Empress® CAD Blocks

Ivoclar Vivadent

Highly aesthetic, leucite-reinforced, glass-ceramic blocks. Biaxial flexural strength of 185 MPa. For fabricating single-unit restorations. Easy-polish option for efficiency or glaze fire for increased strength. Furnace required postmilling.



HT (High Translucency)

Ideal for fabricating smaller restorations, inlays, onlays, and veneers. Provide a lifelike fluorescence and surround the dentition to refine the shade of restoration, achieving a balanced, true-to-nature, chameleon-like effect.

LT (Low Translucency)

Characterized by their high level of brightness and true-to-nature chroma. Suitable for fabricating larger restorations (partial crowns, anterior and posterior crowns, and veneers).

Multi (Multi-chromatic)

Creates a smooth shade transition between dentin and incisal areas to produce highly aesthetic results. Also has a high level of chroma and lifelike opacity in cervical/dentin areas. Ideal for fabricating anterior and posterior crowns, partial crowns, and veneers.

For use with **CEREC®** and **inLab®**



HT Blocks

Size	I10	I12
Shade		
100	—	—
200	—	—
300	—	—
A1	(945-1323)	(945-1335)
A2	(945-1326)	(945-1337)
A3	(945-1327)	(945-1340)
A3.5	(945-1328)	(945-1341)
B1	(945-1329)	(945-1342)
B2	(945-1330)	(945-1343)
B3	(945-1331)	(945-1344)
C2	(945-1332)	(945-1345)
D3	(945-1334)	(945-1346)
5 per Box		

LT Blocks

Size	I12	C14
Shade		
100	—	—
200	—	—
300	—	—
A1	(945-1426)	(945-1437)
A2	(945-1427)	(945-1439)
A3	(945-1429)	(945-1440)
A3.5	(945-1430)	(945-1441)
B1	(945-1431)	(945-1442)
B2	(945-1433)	(945-1443)
B3	(945-1434)	(945-1444)
BL1	(945-1452)	(945-1456)
BL2	(945-1453)	(945-1458)
BL3	(945-1454)	(945-1459)
BL4	(945-1455)	(945-1460)
C2	(945-1435)	(945-1445)
D3	(945-1436)	(945-1446)
5 per Box		

Multi Blocks

Size	I12	C14	C14L
Shade			
A1	(945-1461)	(945-1466)	(945-1478)
A2	(945-1462)	(945-1468)	(945-1483)
A3	(945-1463)	(945-1469)	(945-1485)
A3.5	(945-1464)	(945-1470)	(945-1486)
B1	(945-1465)	(945-1472)	(945-1487)
BL1	(945-1474)	(945-1476)	(945-1488)
BL3	(945-1475)	(945-1477)	(945-1489)
5 per Box		5 per Box	

IPS Empress® CAD Blocks (Cont'd)

For use with Planmeca® PlanMill®



HT Blocks

Size	I10	I12
Shade		
A1	(945-0676)	(945-0688)
A2	(945-0679)	(945-0689)
A3	(945-0680)	(945-0690)
A3.5	(945-0681)	(945-0691)
B1	(945-0682)	(945-0692)
B2	(945-0684)	(945-0693)
B3	(945-0685)	(945-0694)
C2	(945-0686)	(945-0695)
D3	(945-0687)	(945-0696)
	5 per Box	

LT Blocks

Size	I12	C14
Shade		
A1	(945-0710)	(945-0720)
A2	(945-0711)	(945-0721)
A3	(945-0712)	(945-0722)
A3.5	(945-0713)	(945-0723)
B1	(945-0714)	(945-0724)
B2	(945-0715)	(945-0725)
B3	(945-0716)	(945-0726)
C2	(945-0717)	(945-0727)
D3	(945-0719)	(945-0729)
BL1	(945-0736)	(945-0740)
BL2	(945-0737)	(945-0741)
BL3	(945-0738)	(945-0742)
BL4	(945-0739)	(945-0743)
	5 per Box	

Multi Blocks

Size	I12	C14	C14L
Shade			
A1	(945-0879)	(945-0887)	(945-0871)
A2	(945-0880)	(945-0888)	(945-0872)
A3	(945-0882)	(945-0889)	(945-0873)
A3.5	(945-0883)	(945-0891)	(945-0875)
B1	(945-0884)	(945-0893)	(945-0876)
BL1	(945-0885)	(945-0894)	(945-0877)
BL3	(945-0886)	(945-0895)	(945-0878)
	5 per Box		5 per Box

NEW!

For use with PrograMill®



HT Blocks

Size	I12
Shade	
A1	(945-1921)
A2	(945-1925)
A3	(945-1926)
A3.5	(945-1927)
B1	(945-1928)
B2	(945-1930)
B3	(945-1931)
C2	(945-1932)
D3	(945-1933)
	5 per Box

LT Blocks

Size	C14
Shade	
A1	(945-1908)
A2	(945-1911)
A3	(945-1912)
A3.5	(945-1913)
B1	(945-1914)
B2	(945-1915)
B3	(945-1916)
C2	(945-1918)
D3	(945-1920)
BL1	(945-1904)
BL2	(945-1905)
BL3	(945-1906)
BL4	(945-1907)
	5 per Box



IPS e.max® ZirCAD Blocks

Ivoclar Vivadent

Zirconium-oxide ceramic, preshaded, monochromatic blocks. High flexural strength of 850–1,200 MPa. Lifelike aesthetics and fracture toughness. Suitable for single-visit treatments in the dental practice. For full-contour crowns and bridges. No risk of chipping. Furnace required postmilling.



For use with **CEREC®** and **inLab®**



LT Blocks

Size	C17	B45
Shade		
A1	(945-1641)	(945-1648)
A2	(945-1644)	(945-1650)
A3	(945-1645)	(945-1651)
B1	(945-1647)	(945-1652)
B2	(945-1699)	(945-1704)
C2	(945-1700)	(945-1705)
D2	(945-1701)	(945-1706)
Bleach	(945-1698)	(945-1703)
	5 per Box	3 per Box

MT Multi Blocks

Size	C17	B45
Shade		
A1	(945-2083)	(945-2072)
A2	(945-2084)	(945-2074)
A3	(945-2085)	(945-2075)
B1	(945-2087)	(945-2076)
B2	(945-2089)	(945-2079)
C2	(945-2090)	(945-2080)
D2	(945-2091)	(945-2081)
BL1	(945-2082)	(945-2070)
	5 per Box	3 per Box

For use with **Planmeca PlanMill®**



MT Multi Blocks

Size	C17	B45
Shade		
A1	(945-2152)	(945-2161)
A2	(945-2153)	(945-2162)
A3	(945-2154)	(945-2163)
B1	(945-2155)	(945-2164)
B2	(945-2156)	(945-2167)
C2	(945-2157)	(945-2169)
D2	(945-2158)	(945-2170)
BL1	(945-2149)	(945-2160)
	5 per Box	3 per Box



C17

B45

Telio® CAD Temporary Blocks

Ivoclar Vivadent

Acrylate polymer (PMMA) blocks for chairside milling of temporary restorations. Long-term temporaries with a wear period up to 12 months. Flexural strength of 130 MPa. High material homogeneity, durable shade stability, and lifelike fluorescence. Indications include temporary anterior/posterior crown and bridges (with up to 2 pontics), temporary restorations on implants, and therapeutic restorations to correct TMJ problems and occlusal adjustments. Size A16 (S/L) blocks are indicated for temporary, single-tooth, hybrid-abutment crowns and feature a prefabricated interface for direct cementing onto the Sirona® TiBase. No furnace required.



For use with **CEREC®** and **inLab®**

Size	B40L		B55
Shade			
A1	(945-1578)	(945-1588)	(945-1600)
A2	(945-1580)	(945-1589)	(945-1601)
A3	(945-1581)	(945-1590)	(945-1602)
A3.5	(945-1582)	(945-1591)	(945-1605)
B1	(945-1583)	(945-1592)	(945-1606)
BL3	(945-1577)	(945-1586)	(945-1599)
	3 per Box	9 per Box	3 per Box

Abutment Blocks

Size	A16 L	A16 S
Shade		
A1	(945-1615)	(945-1609)
A2	(945-1616)	(945-1610)
A3	(945-1617)	(945-1611)
A3.5	(945-1618)	(945-1612)
B1	(945-1619)	(945-1613)
BL3	(945-1614)	(945-1608)
	3 per Box	



For use with **Planmeca PlanMill®**

Size	B40L		B55
Shade			
A1	(945-0898)	(945-0911)	(945-0925)
A2	(945-0901)	(945-0912)	(945-0926)
A3	(945-0902)	(945-0913)	(945-0927)
A3.5	(945-0905)	(945-0914)	(945-0928)
B1	(945-0907)	(945-0915)	(945-0929)
BL3	(945-0897)	(945-0908)	(945-0924)
	3 per Box	9 per Box	3 per Box



Tetric® CAD

Ivoclar Vivadent

Hybrid ceramic-resin block. High flexural strength of 272 MPa. Restorations blend with the surrounding tooth structure due to the pronounced chameleon-like effect. Indicated for fabricating single-tooth restorations including veneers, inlays, onlays, and crowns. No furnace required: mill, polish, and place.



For use with CEREC® and inLab®

Starter Kit
(945-1985) Ea.

Contains: 9 Tetric CAD blocks (2 of each size I12, HT, shades A1, A2 & A3; & 1 of each size C14 MT, shades A1, A2 & A3), 1–2.5 Gm. syringe of Variolink Esthetic DC, 1–0.5 mL syringe of Adhese Universal VivaPen & 1 OptraPol test pack.

HT Blocks

Size	I12	C14
Shade		
A1	(945-1986)	(945-1997)
A2	(945-1987)	(945-1998)
A3	(945-1988)	(945-1999)
A3.5	(945-1989)	(945-2000)
5 per Box		

MT Blocks

Size	I12	C14
Shade		
A1	(945-1990)	(945-2001)
A2	(945-1992)	(945-2004)
A3	(945-1993)	(945-2005)
A3.5	(945-1995)	(945-2007)
BL	(945-1996)	(945-2009)
5 per Box		

For use with Planmeca PlanMill®

HT Blocks

Size	I12	C14
Shade		
A1	(945-2036)	(945-2050)
A2	(945-2037)	(945-2053)
A3	(945-2038)	(945-2054)
A3.5	(945-2039)	(945-2056)
5 per Box		

MT Blocks

Size	I12	C14
Shade		
A1	(945-2042)	(945-2057)
A2	(945-2043)	(945-2058)
A3	(945-2044)	(945-2059)
A3.5	(945-2045)	(945-2060)
BL	(945-2048)	(945-2061)
5 per Box		



NEW!

For use with Programill®

HT Blocks

Size	I12	C14
Shade		
A1	(945-2294)	(945-2303)
A2	(945-2295)	(945-2305)
A3		(945-2307)
A3.5	(945-2296)	
5 per Box		

MT Blocks

Size	I12	C14
Shade		
A1	(945-2299)	(945-2308)
A2	(945-2300)	(945-2309)
A3	(945-2301)	(945-2310)
5 per Box		



I12 C14

**ivoclar
vivadent**
passion vision innovation



CELTRA® Duo Blocks

Dentsply Sirona

A high-strength, glass-ceramic, CAD/CAM block made of zirconia-reinforced lithium silicate (ZLS). Provides exquisite aesthetics and exceptional strength. Indicated for single-unit posterior and anterior crowns, inlays, onlays, and veneers. Offers 2 processing options: mill and polish (flexural strength of 210 MPa) or mill and fire (flexural strength of 370 MPa).

For use with **CEREC®** and **inLab®**



LT Blocks

Size	C14
Shade	
A1	(222-2164)
A2	(222-2169)
A3	(222-2170)
A3.5	(222-2172)
	4 per Box

Size	C14
Shade	
B1	(222-2173)
B2	(222-2174)
BL2	(222-2175)
BL3	(222-2176)
	4 per Box

Size	C14
Shade	
C1	(222-2177)
C2	(222-2179)
D2	(222-2180)
D3	(222-2181)
	4 per Box

HT Blocks

Size	C14
Shade	
A1	(222-2182)
A2	(222-2183)
A3	(222-2185)
B1	(222-2186)
	4 per Box

Size	C14
Shade	
C1	(222-2187)
C2	(222-2188)
D2	(222-2189)
D3	(222-2190)
	4 per Box

CEREC® Zirconia

Dentsply Sirona

For use with **CEREC®** and **inLab®**



Size	Medi S
Shade	
A1	(587-2365)
A2	(587-2366)
A3	(587-2368)
A3.5	(587-2367)
A4	(587-2369)
B2	(587-2370)
B3	(587-2371)
C2	(587-2372)
C3	(587-2374)
D3	(587-2375)
	Each

Size	Mono L
Shade	
A1	(587-2352)
A2	(587-2353)
A3	(587-2354)
A3.5	(587-2356)
A4	(587-2358)
B2	(587-2359)
B3	(587-2360)
C2	(587-2361)
C3	(587-2362)
D3	(587-2364)
	Each

Size	Meso L
Shade	
A1	(587-8355)
A2	(587-8360)
A3	(587-8356)
A3.5	(587-8359)
	Each

Size	Meso S
Shade	
A1	(587-8358)
A2	(587-8354)
A3	(587-8361)
A3.5	(587-8357)
	Each

inCoris ZI/TZI C

Dentsply Sirona

For use with **CEREC®** and **inLab®**



Size	ZI
Shade	
F0.5 L	(587-2709)
F0.5 S	(587-2708)
F2 L	(587-2711)
F2 S	(587-2710)
20/19 F0.5	(939-5806)
20/19 F0.5	(939-5805)*
	Each

*3 pack

Size	TZI-C Medi
Shade	
A1 3	(587-2400)
A2 3	(587-2401)
A3.5	(587-2403)
A3 3	(587-2402)
A4 3	(587-2405)
B2 3	(587-2406)
B3 3	(587-2407)
C2 3	(587-2408)
C3 3	(587-2409)
D3 3	(587-2435)
	3 Pack

Size	TZI-C Maxi
Shade	
A1 1	(587-2436)
A2 1	(587-2437)
A3 1	(587-2438)
	Each

Shade Guide	
	(587-2351)
	Each

Size	TZI-C Mono
Shade	
A1 3	(587-2384)
A2 3	(587-2387)
A3.5	(587-2389)
A3 3	(587-2388)
A4 3	(587-2390)
B2 3	(587-2392)
B3 3	(587-2393)
C2 3	(587-2396)
C3 3	(587-2398)
D3 3	(587-2399)
	3 Pack



CEREC® Block C/PC

Dentsply Sirona

For use with CEREC® and inLab®

Size	10
Shade	
A1C	(587-2709)
A2C	(587-2708)
A3C	(587-2711)
A3.5C	(587-2710)
A4C	(939-5806)
Bleach 2C	(939-5805)
B2C	(587-2453)
B3C	(587-2454)
	8 Pack

Size	10
Shade	
C2C	(587-2455)
C3C	(587-2456)
D3C	(587-2457)
A1C	(587-2445)
A2C	(587-2446)
Bleach 2C	(587-2458)
	8 Pack

Size	PC-12
Shade	
A1C	(587-2448)
A2C	(587-2449)
A3C	(587-2450)
A3.5C	(587-2451)
	8 Pack

Size	PC-14
Shade	
A1C	(587-2463)
A2C	(587-2465)
A3C	(587-2466)
A3.5C	(587-2467)
A1C	(587-2459)
A2C	(587-2460)
A3C	(587-2461)
A3.5C	(587-2462)
	8 Pack

CEREC® Guide Block

Dentsply Sirona

For use with CEREC® and inLab®

Size	Guide
Maxi	(587-2802)
Medi	(587-2823)
	Each


**Dentsply
Sirona**
VITA SUPRINITY® PC

VITA North America

Partially crystallized, zirconia-reinforced glass ceramic blocks. Biaxial strength of 541 MPa and 3-point flexural strength of 495 MPa. Brilliant esthetics with integrated translucency, opalescence, and fluorescence. Marginal accuracy and high edge stability. Can be milled with very thin edges to prevent chipping. Indicated for posterior and anterior crowns, veneers, and inlays and onlays. Furnace required postmilling.

*For use with CEREC® and inLab®***HT Blocks**

Size	14
Shade	
1M1	(809-0716)
1M2	(809-0717)
2M2	(809-0718)
3M2	(809-0719)
4M2	(809-0720)
	5 per Box

T Blocks


Size	14
Shade	
1M1	(809-0722)
1M2	(809-0723)
2M2	(809-0724)
3M2	(809-0725)
4M2	(809-0726)
	5 per Box


VITA®

VITA ENAMIC® Blocks

VITA North America

Dual-network ceramic material reinforced by a polymer network that delivers resiliency and force-absorption for restorations. Flexural strength of 150–160 MPa. Offers a balanced blend of energy distribution, fracture resistance, and reliability. Excellent bond strength and material properties almost identical to those of natural teeth. Indications include single anterior and posterior full-contour crowns, inlays, onlays, veneers, implant-supported crowns, and implant-abutment restorations. No furnace required: mill, polish, and place.



HT (High Translucency)

Suitable for inlays, onlays, and full or partial crowns. Also used for veneers, cervical veneers, and incisal edges.

T (Translucent)

Suitable for full crowns intended to mask discolored or dark stumps such as in cases of tetracycline discoloration, amalgam tattoos, and metal core buildups.

For use with **CEREC®** and **inLab®**



HT Blocks

Size	EM-10	EM-14
Shade		
0M1		(809-0360)
1M1	(809-0356)	(809-0361)
1M2	(809-0357)	(809-0362)
2M2	(809-0358)	(809-0363)
3M2	(809-0359)	(809-0364)
	5 per Box	5 per Box

T Blocks

Size	EM-14
Shade	
0M1	(809-0366)
1M1	(809-0367)
1M2	(809-0368)
2M2	(809-0369)
3M2	(809-0370)
	5 per Box

For use with **Planmeca PlanMill®**



HT Blocks

Size	EM-14
Shade	
0M1	(809-0219)
1M1	(809-0220)
1M2	(809-0221)
2M2	(809-0222)
3M2	(809-0223)
	5 per Box

VITA ENAMIC® multiColor Blocks

VITA North America

Dual-network ceramic material that features 6 layers of shade gradient, from the translucent incisal edge to the higher chroma layering in the neck. Flexural strength of 150–160 MPa. High-strength, bonded ceramic—combined with the natural properties of polymer—ensure elasticity similar to the natural tooth structure and integrated crack prevention. Indicated for anterior and posterior crowns, implant restorations and veneers. Fast milling with no post-firing required: simply mill, polish, and place.



For use with **CEREC®** and **inLab®**



HT Blocks

Size	EMC-14
Shade	
1M1	(809-0371)
1M2	(809-0372)
2M2	(809-0373)
3M2	(809-0374)
4M2	(809-0375)
	5 per Box

VITA®



VITA CAD/CAM BLOCKS

With over 30 years of CAD/CAM excellence, only VITA offers the widest range of material to perform virtually any restoration procedure chairside.



**BUY MORE,
GET MORE!**

AVAILABLE IN VITA CLASSICAL A1-D4® AND VITA SYSTEM 3D-MASTER®



VITA SUPRINITY® PC

Partially crystallized, zirconia-reinforced block that is strong without compromising esthetics



VITA ENAMIC® multiColor

Reinforced ceramic/polymer dual structure with no firing required, that mimics natural teeth



VITABLOCS® MARK II

High quality feldspathic porcelain — the CAD/CAM gold standard




VITA ENAMIC® IS

Hybrid ceramic for superstructures with high load-bearing capacity



VITABLOCS® TriLuxe forte

Multi-shade layers for highly esthetic anterior restorations



VITA ZIRCONIA YZ® HT

Up to 1200 MPa strength with excellent light-transmitting properties for esthetics

BUY 4 CEREC BLOCKS, GET 1 BLOCKS FREE!*

BUY 12 CEREC BLOCKS, GET 4 BLOCKS FREE!*

BUY 25 CEREC BLOCKS, GET VITA EASYSHADE® V FREE!

BUY 25 CEREC BLOCKS, GET 10 BLOCKS FREE!*

BUY 30 CEREC BLOCKS, GET VITA AKZENT PLUS KIT, CERAMIC ETCH & ENAMIC TECHNICAL POLISHING KIT FREE!

BUY 50 CEREC BLOCKS, GET 20 BLOCKS FREE!*

BUY 100 CEREC BLOCKS, GET 40 BLOCKS FREE!*

*equal or lesser value

800-828-3839
vitanorthamerica.com

VITA – perfect match

VITA


PROMO CODE: Q22021CA. Valid on purchases made April 1-June 30, 2021. Proof of purchase on one invoice required for free goods. Invoice must be submitted to promo@vitanorthamerica.com no later than July 15, 2021. Offers valid in Canada only. Offers may not be combined with any other promotion or discounts. VITA North America reserves the right to discontinue or change the offers at any time without notice. No cash value.

© 2021 VITA North America. VITA SUPRINITY® PC, VITA ENAMIC®, VITABLOCS® and VITA ZIRCONIA YZ® are registered trademarks of VITA Zahnfabrik

VITA ENAMIC® IS Implant Solutions CAD/CAM Blocks

VITA North America

Dual-network ceramic material with all of the characteristics of a natural tooth. Flexural strength of 150–160 MPa. Delivers unrivaled resiliency and force-absorption for long-term implant success. Provides a balanced blend of energy distribution, fracture toughness, and reliability. Indicated for abutment crowns (IS-16) and milling abutments (IS-14). Integrated interface for Sirona® TiBase. No furnace required.



For use with **CEREC®** and **inLab®**

HT Blocks

Size	IS-16L	IS-16S
Shade		
1M1	(809-0517)	(809-0522)
1M2	(809-0518)	(809-0523)
2M2	(809-0519)	(809-0524)
3M2	(809-0520)	(809-0525)
4M2	(809-0521)	(809-0526)
	5 per Box	

T Blocks

Size	IS-14L	IS-14S
Shade		
1M1	(809-0527)	(809-0532)
1M2	(809-0528)	(809-0533)
2M2	(809-0529)	(809-0534)
3M2	(809-0530)	(809-0535)
4M2	(809-0531)	(809-0536)
	5 per Box	

VITABLOCS® Mark II Blocks

VITA NORTH AMERICA

Fine-structure, monochromatic, feldspar-ceramic material. Flexural strength of 154 (±15) MPa. Excellent translucency enables a natural play of colors. Mimics natural aesthetics. Indications include anterior/posterior crowns, inlays, onlays, and veneers. Reliable and easy: no furnace required, just mill and polish.



For use with **CEREC®** and **inLab®**

Classical Shades

Size	I10	I12	I14
Shade			
A1C	(809-0379)	(809-0388)	(809-0399)
A2C	(809-0380)	(809-0389)	(809-0400)
A3C	(809-0381)	(809-0390)	(809-0401)
A3.5	(809-0382)	(809-0391)	(809-0402)
A4	(809-0383)	(809-0392)	(809-0403)
B2	(809-0384)	(809-0393)	(809-0404)
B3	(809-0385)	(809-0394)	(809-0405)
C2	(809-0386)	(809-0395)	(809-0406)
C3	(809-0398)	(809-0396)	(809-0407)
D3	(809-0387)	(809-0397)	(809-0408)
	5 per Box		5 per Box

3D Master Shades

Size	I10	I12	I14
Shade			
0M1C	-	(809-0493)	(809-0409)
1M1C	(809-0484)	(809-0494)	(809-0410)
1M2C	(809-0485)	(809-0495)	(809-0411)
2M1C	(809-0486)	(809-0496)	(809-0412)
2M2C	(809-0487)	(809-0497)	(809-0413)
2M3C	(809-0488)	(809-0498)	(809-0414)
3M1C	(809-0489)	(809-0499)	(809-0415)
3M2C	(809-0490)	(809-0500)	(809-0416)
3M3C	(809-0491)	(809-0501)	(809-0417)
4M2C	(809-0492)	(809-0502)	(809-0418)
	5 per Box		5 per Box

VITA®



VITABLOCS® TriLuxe forte Blocks

VITA North America

Multichromatic feldspar-ceramic material. Flexural strength of 154 (± 15) MPa. Mimics natural aesthetics. Features 4 gradient layers with smooth transitions from neck to incisal edge. Indicated for anterior crowns and veneers. No furnace required: just mill, polish, and place.



For use with **CEREC®** and **inLab®**



Size	TF-12	TF-14/12	TF-14/14
Shade			
1M2C	(809-0431)	(809-0438)	(809-0537)
2M2C	(809-0432)	(809-0439)	(809-0538)
3M2C	(809-0433)	(809-0440)	(809-0539)
A1	(809-0434)	(809-0441)	(809-0540)
A2	(809-0435)	(809-0442)	(809-0541)
A3	(809-0436)	(809-0443)	(809-0542)
A 3.5	(809-0437)	(809-0444)	(809-0543)
	5 per Box	5 per Box	5 per Box

For use with **Planmeca PlanMill®**



Size	TF-14/14
Shade	
A1C	(809-0242)
A2C	(809-0243)
A3C	(809-0241)
A3.5C	(809-0244)
	5 per Box

VITABLOCS® RealLife Blocks

VITA North America

Fine-structure, multichromatic, feldspar-ceramic material. Flexural strength of 154 (± 15) MPa. Offers 3-dimensional aesthetics utilizing the Sirona® software solution. Creates detailed simulation of natural tooth structure and aesthetics. Replicates the curved color gradient found in natural dentition. Eliminates the need for characterization with stains or individualizing with layered veneering materials. Indications include full and partial anterior crowns and veneers, up to and including premolars. No furnace required.




For use with **CEREC®** and **inLab®**



HT Blocks

Size	RL-14/14
Shade	
0M1	(809-0555)
1M1	(809-0556)
1M2	(809-0552)
2M1	(809-0557)
2M2	(809-0553)
3M2	(809-0554)
	5 per Box



VITA®

VITA YZ® HT Blocks

VITA North America

Highly translucent zirconia material. High strength (approx. 1,200 MPa) and excellent light-refraction properties.

High edge stability and precise sintering results ensure accurate restorations. Reduction of wall thickness enables thinner restorations and preserves tooth structure.

Indications include monolithic anterior and posterior bridge restorations and in cases of limited available space. Noncolored (white) blocks for individual coloring with VITA YZ® HT shade liquids. Furnace required postmilling.

For use with **CEREC®** and **inLab®**



HT Blocks

Size	20/19	40/19	55/19	65/25
Shade				
White	(809-0510)	(809-0514)	(809-0515)	(809-0516)
	4 per Box	2 per Box	2 per Box	Each



VITA CAD-TEMP™ monoColor Blocks

VITA North America

Cross-linked, monochromatic acrylate-polymer material that features a fiber-free, homogeneous, and high-molecular structure with microfiller. Flexural strength of >80 MPa. Well-balanced combination of tensile strength and elasticity. Indicated for fabricating multiunit, fully or partially anatomical, long-term temporary bridge restorations (up to 2 pontics). Bridge block is also suitable for temporary crown restorations. Clinical wearing period of at least 2 years. No furnace required.

For use with **CEREC®** and **inLab®**



Size	CT-40		CT-55
Shade			
0M1T	(809-0506)	(809-0558)	—
1M2T	(809-0507)	(809-0559)	(809-0448)
2M2T	(809-0508)	(809-0560)	(809-0449)
3M2T	(809-0509)	(809-0562)	(809-0450)
	2 per Box	10 per Box	Each



VITA CAD-TEMP™ multiColor Blocks

VITA North America


Cross-linked, acrylate-polymer material with 4 layers of varying shade intensity. Features a fiber-free, homogeneous, and high-molecular structure with microfiller. Flexural strength of >80 MPa. Well-balanced combination of tensile strength and elasticity. Indicated for fabricating multiunit, fully or partially anatomical, long-term temporary bridge restorations (up to 2 pontics). Bridge block is also suitable for temporary crown restorations. Clinical wearing period up to 3 years. No furnace required.

For use with **CEREC®** and **inLab®**



Size	CTM-40	
Shade		
1M2T	(809-0511)	(809-0445)
2M2T	(809-0512)	(809-0446)
3M2T	(809-0513)	(809-0447)
	2 per Box	10 per Box

VITA®



KATANA™ Zirconia Blocks

Kuraray America

Super-translucent, multilayered (STML) zirconia block with mechanical properties and aesthetic qualities that mimic natural teeth. Flexural strength of 763 MPa. Dry milling is recommended. Allows a 18-min. sinter time with CEREC® Speedfire. Indicated for fabricating ceramic restorations: anterior and posterior crowns, inlays, onlays, and veneers. Furnace required postmilling.



For use with **CEREC®** and **inLab®**

HT Blocks

Introductory Sets

Natural

(721-0882) Ea.
Contains: 2-12Z blocks (1 of each shade: A2 & A3).

White

(721-0883) Ea.
Contains: 2-12Z blocks (1 of each shade: NW & B1).

Basic Kit

(721-0881) Ea.
Contains: 5-12Z blocks (1 of each shade: A1, A2, A3, A3.5 & B2).

Size	12Z	12Z	14Z	14Z	14Z L	14Z L
Shade						
A1	(721-0851)	-	(721-0858)	-	(721-1080)	-
A2	(721-0852)	-	(721-0859)	-	(721-1081)	-
A3	(721-0853)	-	(721-0860)	-	(721-1082)	-
A3.5	(721-0854)	-	(721-0861)	-	(721-1083)	-
A4	(721-0855)	-	(721-0862)	-	-	(721-1085)
B1	(721-0856)	-	(721-0863)	-	(721-1084)	-
B2	(721-0857)	-	(721-0864)	-	-	(721-1086)
NW	(721-0879)	-	(721-0880)	-	-	(721-1093)
B3	-	(721-0865)	-	(721-0871)	-	(721-1087)
C1	-	(721-0866)	-	(721-0872)	-	(721-1088)
C2	-	(721-0867)	-	(721-0874)	-	(721-1089)
C3	-	(721-0868)	-	(721-0875)	-	(721-1090)
D2	-	(721-0869)	-	(721-0876)	-	(721-1091)
D3	-	(721-0870)	-	(721-0877)	-	(721-1092)
CL*	-	(721-0850)	-	-	-	-
	5 per Box	2 per Box	5 per Box	2 per Box	3 per Box	2 per Box

*Not a multilayered block

KATANA™ AVENCIA™ Blocks

Kuraray America

Hybrid ceramic composite resin material. Flexural strength of 230 MPa to withstand occlusal pressure. Excellent wear resistance, polishability, and gloss retention. Indicated for the fabrication of inlays, onlays, veneers, and full-contour crowns. No furnace required: mill, polish, and place.



For use with **CEREC®** and **inLab®**

Size	12	14L
Shade		
A2	(721-1032)	(721-1035)
A3	(721-1033)	(721-1036)
A3.5	(721-1034)	(721-1037)
	5 per Box	5 per Box



kuraray

CERASMART™ Blocks

GC America

Force-absorbing, flexible nanoceramic CAD/CAM material. High, 3-point flexural strength of 231 MPa and a biaxial strength of 238 MPa. Combines the best characteristics of a high-strength ceramic and a composite. Provides the highest degree of flexibility, strength, and breaking energy, as well as superior reliability, aesthetics, milling accuracy, and self-polishing properties. Additionally, it offers properties almost identical to those of natural teeth. Can be used for posterior, anterior, inlay, onlay, and implant restorations. Also enables only the slightest amount of tooth reduction for minimally invasive restorations. No furnace required.

**HT (High Translucency)**

For replacing enamel. Use in cases of inlays, onlays, veneers, partials, and full crowns.

LT (Low Translucency)

For replacing dentin and enamel structures. Used to mask discolored preparations, particulary crowns.

For use with **CEREC®** and **inLab®**

**HT Blocks**

Size	12	14	14L
Shade			
A1	(333-4437)	(333-4493)	(333-4549)
A2	(333-4442)	(333-4496)	(333-4550)
A3	(333-4445)	(333-4497)	(333-4551)
A3.5	(333-4447)	(333-4498)	(333-4552)
B1	(333-4448)	(333-4508)	(333-4553)
	5 per Box	5 per Box	5 per Box

LT Blocks

Size	12	14	14L
Shade			
A1	(333-4450)	(333-4509)	(333-4554)
A2	(333-4452)	(333-4539)	(333-4555)
A3	(333-4453)	(333-4540)	(333-4556)
A3.5	(333-4454)	(333-4541)	(333-4557)
B1	(333-4455)	(333-4542)	(333-4559)
	5 per Box	5 per Box	5 per Box

Bleach Shades

Size	12	14	14L
Shade			
Bleach	(333-4491)	(333-4548)	(333-4560)
	5 per Box	5 per Box	5 per Box

For use with **Planmeca PlanMill®**

**HT Blocks**

Size	12	14	14L
Shade			
A1	(333-4127)	(333-4167)	(333-4144)
A2	(333-4128)	(333-4169)	(333-4145)
A3	(333-4129)	(333-4170)	(333-4146)
A3.5	(333-4130)	(333-4171)	(333-4147)
B1	(333-4131)	(333-4172)	(333-4149)
	5 per Box	5 per Box	5 per Box

LT Blocks

Size	12	14	14L
Shade			
A1	(333-4132)	(333-4173)	(333-4150)
A2	(333-4133)	(333-4174)	(333-4151)
A3	(333-4136)	(333-4175)	(333-4154)
A3.5	(333-4139)	(333-4176)	(333-4162)
B1	(333-4140)	(333-4178)	(333-4165)
	5 per Box	5 per Box	5 per Box

Bleach Shades

Size	12	14	14L
Shade			
Bleach	(333-4141)	(333-4179)	(333-4166)
	5 per Box	5 per Box	5 per Box



12

14

14L

CERASMART® 270 **New!**

GC America

Esthetic hybrid block made to provide fast and long-lasting indirect restorations. The newly developed Full-Coverage Silane Coating (FSC) with improved nanofiller technology offers high physical and esthetic properties. Milling is fast and efficient, making it possible to create a qualitative restoration with fine margins and optimal fit in a very short time. It may also be polished or characterized for more vibrant esthetics. Mills quickly and is gentle on milling burs, virtually no chipping and no need for sintering/crystallization firing. Flexible, tough, and fracture resistant. Easy to polish, retains gloss over time, and self polishes.



For use with **CEREC®** and **inLab®**

HT Blocks

Size	12	14	14L
Shade			
A1	(333-5258)	(333-5272)	(333-5246)
A2	(333-5260)	(333-5275)	(333-5247)
A3	(333-5261)	(333-5277)	(333-5248)
A3.5	(333-5262)	(333-5278)	(333-5249)
B1	(333-5264)	(333-5280)	(333-5250)
	5 per Box	5 per Box	5 per Box

LT Blocks

Size	12	14	14L
Shade			
A1	(333-5265)	(333-5281)	(333-5251)
A2	(333-5267)	(333-5284)	(333-5252)
A3	(333-5268)	(333-5285)	(333-5254)
A3.5	(333-5269)	(333-5286)	(333-5255)
B1	(333-5270)	(333-5287)	(333-5256)
	5 per Box	5 per Box	5 per Box

Bleach

Size	12	14	14L
Shade			
Bleach	(333-5271)	(333-5288)	(333-5257)
	5 per Box	5 per Box	5 per Box



GC Initial™ LRF Blocks

GC America

Leucite-reinforced feldspar blocks for full-contour, all-ceramic indirect restorations. Flexural strength of 210 MPa when polished and 250 MPa after using glaze paste. High fracture stability. Natural, chameleon-like effect. Indicated for single anterior and posterior partial and full crowns, inlays and onlays, veneers, and endo crowns on molars. Does not require firing: just mill, polish, and cement.



For use with **CEREC®** and **inLab®**

HT Blocks


Size	12	14	14L
Shade			
A1	(333-3976)	(333-4060)	(333-4570)
A2	(333-4577)	(333-4579)	(333-4580)
A3	(333-4586)	(333-4589)	(333-4595)
A3.5	(333-4604)	(333-4605)	(333-4606)
B1	(333-4612)	(333-4613)	(333-4614)
	5 per Box		

LT Blocks

Size	12	14	14L
Shade			
A1	(333-4572)	(333-4573)	(333-4574)
A2	(333-4582)	(333-4583)	(333-4585)
A3	(333-4597)	(333-4600)	(333-4602)
A3.5	(333-4607)	(333-4610)	(333-4611)
B1	(333-4615)	(333-4616)	(333-4618)
	5 per Box		

Bleach Blocks

Size	12	14	14L
Shade			
Bleach	(333-4619)	(333-4620)	(333-4621)
	5 per Box		



GC

3M™ Chairside Zirconia

3M

Preshaded zirconia blocks. High flexural strength of >800 MPa. Low minimum wall thickness of just 0.8mm to preserve more tooth structure. Highly esthetic for natural-looking restorations. Fast sintering time of approximately 20 minutes for a thin-walled crown using CEREC® Speedfire. Indicated for single-unit crowns and 3-unit bridges. Furnace required postmilling.



For use with **CEREC®** and **inLab®**

Size	20	39
Shade		
A1	(777-2434)	(777-2443)
A2	(777-2435)	(777-2444)
A3	(777-2436)	(777-2445)
A3.5	(777-2437)	(777-2446)
B1	(777-2438)	(777-2447)
C1	(777-2440)	(777-2448)
D2	(777-2441)	(777-2449)
BL	(777-2433)	(777-2442)
	5 per Box	3 per Box

3M Science.
Applied to Life.™



RELY ON US FOR ON-SITE EQUIPMENT AND TECHNICAL SERVICE

- Equipment installation and relocation
- Manufacturer warranty support
- Hi-tech support & services
- Prompt and reliable repair services
- Preventative maintenance agreements

1-888-818-5838

 **HENRY SCHEIN®** 
TECHNICAL SERVICES

Rely on Us™

Lava™ Ultimate CAD/CAM Restorative Blocks

3M

Resin-nanoceramic, restorative CAD/CAM material. Flexural strength of 200 MPa. Indicated for permanent adhesive, single-unit restorations including inlays, onlays, and veneers. Fast and easy to mill. Brilliant aesthetics with long-lasting polish. Less wear to opposing enamel and less brittle than glass ceramics. No firing necessary.



For use with **CEREC®** and **inLab®**

Introductory Kit

(777-2342) Ea.

Contains: 10 blocks total (5 size 12 blocks, 1 of each shade: A1/LT, A2/HT, A2/LT, A3/HT & A3/LT; & 5 size 14L blocks, 1 of each shade: A1/LT, A2/HT, A2/LT, A3/HT & A3/LT).

HT Blocks

Size	12	14L
Shade		
A1	(777-2355)	(777-2380)
A2	(777-2358)	(777-2382)
A3	(777-2364)	(777-2384)
B1	(777-2372)	(777-2387)
5 per Box		

LT Blocks

Size	12	14L
Shade		
A1	(777-2356)	(777-2381)
A2	(777-2361)	(777-2383)
A3	(777-2366)	(777-2385)
A3.5	(777-2371)	(777-2386)
B1	(777-2374)	(777-2388)
C2	(777-2375)	(777-2389)
D2	(777-2376)	(777-2390)
BL	(777-2354)	(777-2377)
5 per Box		

For use with **Planmeca PlanMill®**



HT Blocks

Size	14L
Shade	
A1	(777-2196)
A2	(777-2198)
A3	(777-2200)
B1	(777-2205)
5 per Box	

LT Blocks

Size	14L
Shade	
A1	(777-2197)
A2	(777-2199)
A3	(777-2201)
A3.5	(777-2202)
B1	(777-2206)
BL	(777-2195)
C2	(777-2207)
D2	(777-2208)
5 per Box	



12

14L

3M Science.
Applied to Life.™

Dealing with office technology can be
FRUSTRATING!
We can help!

1-800-737-5029 ext 3795 • hsdssupport@henryschein.ca



HENRY SCHEIN®
DIGITAL SOLUTIONS



Henry Schein Cares

Helping Health Happen®

MAKING A DIFFERENCE

Together, we are Helping Health Happen®



Back to School

Providing hundreds of backpacks filled with school supplies and clothing annually to give underprivileged children a sense of confidence as they head back to school.



Community

Helping underserved people in communities across Canada receive dental care they normally could not afford.



Disaster Relief

Shipping more than \$1 million of healthcare products a year to aid when natural disasters strike globally.



Going Green

Reducing paper, planting trees and saving energy. Our Eastern Distribution Centre has gone paperless, eliminating 17,000 lbs of paper usage each year!



Henry Schein Cares

Helping Health Happen®



1-800-668-5558 | www.henryschein.ca



Grandio® blocs

VOCO

Resin nanoceramic hybrid blocks. Biaxial flexural strength of 333 MPa. Ideal aesthetics with restorations that resemble natural teeth perfectly. High filler content for better wear resilience and long-term color stability. Allows for precision milled, thin, tapered edges without the concern for chipping and breakage. Indications include crowns, inlays/onlays, veneers, and implant-supported crowns. No firing required. Complete Kit

(999-2329).....Ea.

Contains: 2 size 12 blocks (1 of each: A2/LT & A3/HT), 3 size 14L blocks (1 of each: A3/LT, A3/HT & A3.5/LT), 1–10 Gm. Bifix QM QuickMix syringe (universal shade), 5–0.1 mL single doses of Futurabond® U, 1–5 mL bottle of Ceramic Bond, 1 Dimanto® complete kit, 10 Single Tim application brushes, 5 mixing tips (type 15) & 5 intraoral tips (type 1).

Assorted Packages

Size 12
(999-2328)..... 5 per Pkg.

Contains: 5 assorted size 12 blocks (1 of each LT shade: A1, B1, C2 & BL; & 1 HT shade A1).

Size 14L
(999-2314) 5 per Pkg.

Contains: 5 assorted size 14L blocks (1 of each LT shade: A1, B1, C2 & BL; & 1 HT shade A1).



For use with CEREC® and inLab®

HT Blocks

Size	12	14L
Shade		
A1	(999-2316)	(999-2295)
A2	(999-2318)	(999-2298)
A3	(999-2320)	(999-2302)
A3.5	(999-2323)	(999-2304)
	5 per Box	5 per Box

LT Blocks

Size	12	14
Shade		
A1	(999-2317)	(999-2297)
A2	(999-2319)	(999-2299)
A3	(999-2321)	(999-2300)
A3.5	(999-2322)	(999-2306)
B1	(999-2325)	(999-2308)
C2	(999-2327)	(999-2313)
BL	(999-2326)	(999-2310)
	5 per Box	5 per Box



VOCO
THE DENTALISTS

Grandio® blocs are not included in CEREC® Software.

**It's Quick, It's Easy, It's Accurate**

With our Aruba eZ™ software system you can build your own tooth orders throughout the workday and upload the completed order to the Zahn Canada® Tooth Department for fast delivery.



ELECTRONIC
TOOTH
ORDERING



ZAHN CANADA®
A DIVISION OF HENRY SCHEIN®

1-800-496-9500  www.zahncanada.ca



DS Universal Stains & Overglazes

Dentsply Sirona

Universal stain and glaze system for customizing restorations including CELTRA® Duo and CEREC® Zirconia.

Stain Refill..... 5 Gm.

Specify:

Body Stains

S0	(811-0702)
S1	(811-0703)
S2	(811-0704)
S3	(811-0705)
S4	(811-0708)

Incisal Stains

Incisal 1	(811-0717)
Incisal 2	(811-0718)

Stains

White	(811-0719)
Crème	(811-0720)
Sunset	(811-0727)
Copper	(811-0728)
Khaki	(811-0729)
Olive	(811-0730)
Mahogany	(811-0731)
Violet	(811-0732)
Raspberry	(811-0733)
Purple	(811-0734)
Grey	(811-0736)
Chestnut	(811-0738)
Blue	(811-0739)
Pink	(811-0740)

Overglaze Refill..... 5 Gm.

Specify:

Overglaze	(811-0741)
Overglaze-High Fluorescent	(811-0742)



IPS e.max® CAD Shade Guide

Ivoclar Vivadent

Shade Guide.....Ea.

Specify:

MT Blocks	(945-0930)
Impulse Blocks	(945-0418)



OPTIGLAZE™ Color

GC America

Light-cure characterization coating developed specifically to stain and glaze CERASMART™ milled blocks, composite, acrylic, and PMMA-type materials.

Complete Kit

(333-3725)Ea.

Contains: 1 of each: 2.6 mL bottle of each shade (A Plus, B Plus, C Plus, Blue, Clear, Clear HV, Grey, Ivory White, Lavender, Olive, Orange, Pink, Pink Orange, Red, Red Brown, White & Yellow), 1 disposable palette, 1 fiat brush, 1 round brush & 1 brush holder.

Color Shades-2.6 mL..... Ea.

Specify:

A Plus	(333-3724)
B Plus	(333-3734)
C Plus	(333-3735)
Blue	(333-3731)
Grey	(333-3743)
Ivory White	(333-3733)
Lavender	(333-3746)
Olive	(333-3742)
Orange	(333-3741)
Pink	(333-3719)
Pink Orange	(333-3722)
Red	(333-3720)
Red Brown	(333-3723)
White	(333-3737)
Yellow	(333-3740)

Color Shades-5 mL..... Ea.

Specify:

Clear	(333-3748)
Clear HV	(333-3749)

OPTIGLAZE™-15 mL

Light-Cure

(333-3259)Ea.



IPS e.max CAD Characterization Kit

Ivoclar Vivadent

Shades, stains, and glazes for enhancing IPS e.max CAD monolithic crowns, abutment crowns, inlays/onlays, and bridges.

(945-1667)Ea.

Contains: 1-3 Gm. syringe of IPS e.max CAD Crystall./Shades (1 of each shade: 1, 2 & 3), 1-3 Gm. syringe of IPS e.max CAD Crystall./Stains (1 of each stain: white, sunset, khaki & mahogany), 1-3 Gm. syringe of IPS e.max CAD Crystall./Glaze Fluo, 1-270 mL bottle of IPS e.max CAD Crystall./Glaze Spray, 1-15 mL bottle of IPS e.max CAD Stain & Glaze Liquid, 1-5 Gm. bottle of IPS e.max CAD Add-on Powder, 1-15 mL bottle of IPS e.max CAD add-on liquid, 1-15 mL bottle of IPS e.max CAD Crystall./Add-On Liquid Allround, 1 crystallization tray, 1 speed crystallization tray & 3 crystallization pins (small, medium & large).



IPS e.max CAD Crystall Shades & Stains

Ivoclar Vivadent

Used to characterize e.max IPS CAD restorations. Can be applied with glaze to the precrystallized restoration to allow for a single firing cycle.

Shades Syringe Refill..... 3 Gm.

Specify:

0	(945-1929)	4	(945-0173)
1	(945-4763)	11	(945-2748)
2	(945-7600)	12	(945-6852)
3	(945-3275)		

Stains Syringe Refill..... 1 Gm.

Specify:

White	(945-0732)
Crème	(945-3975)
Sunset	(945-8690)
Copper	(945-4279)
Olive	(945-0303)
Khaki	(945-9953)
Mahogany	(945-9602)

Accessories

Glaze Paste Syringe

(945-6632) 3 Gm.

(945-1669) 3 Gm.

Glaze Spray

(945-9726)270 mL

Stain and Glaze Liquid

Allows you to optimize the consistency of IPS e.max CAD Crystall stains, shades, and glazes prior to application.

(945-5331) 15 mL

Add-On Powder 5 Gm.

Specify:

Dentin	(945-0984)
Connect	(945-1653)

Incisal Add-on

(945-1174)5 Gm.

Add-On Liquid

(945-6972)15 mL

IPS Object Fix

Convenient substitute for firing pins or for use with veneer cases. Allows the restoration to be supported on the firing tray in a foam gel that promotes proper air flow around the restoration during the firing process.

IPS Object Fix..... 10 mL

Specify:

Putty	(945-1668)
Flow	(945-8518)



IPS Natural Die Material

Ivoclar Vivadent

Light-cured material that simulates the shade of the preparation. Provides the optimum basis on which shades of the original tooth can be reproduced.

Refill Syringe..... 8 Gm.

Specify:

ND2 (945-1801) ND7 (945-4370)
ND3 (945-5206) ND8 (945-0919)
ND4 (945-3982)
ND5 (945-0443)

(945-0102)80 per Pkg.

Separator

(945-7703)20 mL

(623-0065) Ea.



IPS Empress® CAD Add-On and Build-Up

Ivoclar Vivadent

Allows for the addition of contacts on milled IPS Empress CAD restorations.

Add-On Material

(945-7202)20 Gm.



Build-Up Liquid

(945-3667)60 mL

Enamelite™ Ceramic Spray Glaze

Keystone Industries

Stain and glaze in one firing rather than firing a ceramic glaze multiple times. Guarantees shade-matching under natural and artificial light for zirconia, lithium-disilicate, and lithium-silicate restorations. Spray a thin, uniform, and accurate coat in seconds. Provides faster drying than brush-on applications and will not puddle at margins. Compatible with all dental ceramics.



Low-Fusing

Small Can2 oz.
(206-0030)

Medium Can3.9 oz.
(206-0005)

Low-Fusing Fluorescent

Fluorescence is built into the glaze, eliminating negative results in certain lighting conditions.

Small Can2 oz.
(242-0887)



Cerabien™ ZR FC Paste Stain

Kuraray America

Paste-type stain (surface stain) that allows easy characterization of full zirconia as well as layered and pressed zirconia. Can be baked at 750°C during the fabrication of full-zirconia restorations. Lower temperature than used for conventional external stains.

Basic Kit

(366-7509)Ea.

Contains: 1-5 Gm. jar of clear glaze, 7-3 Gm. jars of paste stain (1 of each shade: dark gray, grayish blue, cervical 2, earth brown, salmon pink, A+ & white), 1-10 mL external stain liquid & 1 color stain guide.

Paste Stain Refills

Translucent

5 Gm. Jar Ea.

Specify:

Glaze (366-7507)
Clear Glaze (366-7508)

Enamel

3 Gm. Jar Ea.

Specify:

Light Gray (366-7479)
Dark Gray (366-7481)
Grayish Blue (366-7483)
Blue (366-7484)

Dentin

3 Gm. Jar Ea.

Specify:

Cervical 1 (366-7497) B+ (366-7468)
Cervical 2 (366-7452) C+ (366-7470)
Cervical 3 (366-7453) D+ (366-7472)
A+ (366-7462)

Enhancement

3 Gm. Jar Ea.

Specify:

White (366-7474)
Mamelon Orange 1 (366-7493)
Mamelon Orange 2 (366-7494)
Earth Brown (366-7454)
Reddish Brown (366-7455)

Adjustment

3 Gm. Jar Ea.

Specify:

Value (366-7475)
Black (366-7478)
Green 1 (366-7488)
Green 2 (366-7489)
Yellow (366-7492)
Salmon Pink (366-7459)
Pink (366-7460)
Red (366-7461)
Fluoro (366-7477)

ES

Great for mixing External Stain and Glazing Powder. Keeps Stain moist to make staining work easier.

(366-0152)Ea.



VITA AKZENT® Plus Paste

VITA North America

Ideal for fast and simple characterization of restorations made from VITABLOCS®. Ready-to-use paste with uniform consistency and homogeneous pigmentation for fast application. Can be mixed with one another to achieve individual shade effects, diluted, or mixed again.

Complete Paste Kit

(809-0247) Ea.

Contains: 14-4 Gm. bottles of Effect Stain Paste, 5-4 Gm. bottles of Body Stain Paste, 1-4 Gm. bottle of Glaze Paste, 1-4 Gm. bottle of Glaze LT Paste, 1-4 Gm. bottle of Finishing Agent Paste & 1-15 mL bottle of Paste Fluid.

Effect Stain Paste 4 Gm.

Specify:

White (809-0266)
Lemon (809-0271)
Sunshine (809-0274)
Orange (809-0276)
Russet (809-0278)
Khaki (809-0280)
Pink (809-0282)
Dark Red (809-0284)
Lilac (809-0286)
Blue (809-0287)
Grey-Blue (809-0289)
Grey (809-0291)
Black (809-0293)

Body Stain Paste 4 Gm.

Specify:

Yellow (809-0251)
Yellow-Brown (809-0254)
Orange (809-0257)
Olive-Grey (809-0260)
Grey-Brown (809-0263)

Glaze Paste

(809-0191) 4 Gm.

Glaze LT Paste

(809-0295) 4 Gm.

Finishing Agent Paste

(809-0296) 4 Gm.

Paste Fluid

(809-0249)15 mL



VITA® Classical Shade Guide

VITA North America

Provides an easy-to-use, precise method of matching the various shades of vital, living teeth. Autoclavable.

Shade Guide with Shade Tabs

(809-0106)Ea.

Shade Tab Refill..... Ea.

Specify:

A1 (809-0110) C1 (809-0120)
A2 (809-0111) C2 (809-0121)
A3 (809-0112) C3 (809-0123)
A3.5 (809-0113) C4 (809-0124)
A4 (809-0114) D2 (809-0125)
B1 (809-0115) D3 (809-0126)
B2 (809-0116) D4 (809-0127)
B4 (809-0118)



VITABLOCS® 3D-Master® Shade Guide

VITA North America

Shade tabs fabricated from the same ceramic materials used in the original VITABLOCS. Improves shade selection by matching shades faster and more precisely.

VITABLOCS® 3D-Master® Shade Guide
(809-0245)Ea.



VITA AKZENT® Plus Powder

VITA North America

Stain that features good stability characteristics and shade stability. Ideal for characterizing conventional veneering ceramics as well as the fine-structured, feldspar-ceramic VITA PM 9. Can be mixed freely with each other to obtain more custom shade effects.

Complete Powder Kit
(809-0246)Ea.

Contains: 14-3 Gm. bottles of Effect Stain Powder, 5-3 Gm. bottles of Body Stain Powder, 4-3 Gm. bottles of Chroma Stain, 1-5 Gm. bottle of Glaze, 1-5 Gm. bottle of Glaze LT, 1-5 Gm. bottle of Finishing Agent Powder & 1-20 mL bottle of powder fluid.

Effect Stain Powder 3 Gm.

Specify:
Cream (809-0269)
Lemon (809-0270)
Sunshine (809-0273)
Orange (809-0275)
Russet (809-0277)
Khaki (809-0279)
Pink (809-0281)
Dark Red (809-0283)
Lilac (809-0285)
Grey-Blue (809-0288)
Grey (809-0290)
Black (809-0292)
Blue (809-0186)
White (809-0187)

Body Stain Powder 3 Gm.

Specify:
Yellow (809-0250)
Yellow-Brown (809-0253)
Orange (809-0256)
Olive-Grey (809-0259)
Grey-Brown (809-0262)

Glaze Powder
(809-0160) 5 Gm.

Glaze Powder
(809-0294) 30 Gm.

Glaze LT Powder
(809-0188) 5 Gm.

Powder Fluid
(809-0182) 20 mL

Finishing Agent Powder
(809-0297) 5 Gm.



VITA AKZENT® Plus Spray

VITA North America

Ideal spray for staining large surfaces, especially for monolithic restorations. Ready to use and easy to apply. Spray head features an even application for an easy, uniform surface treatment and a precise spray application without any waste from scatter.

Complete Spray Kit
(809-0248)Ea.

Contains: 5-75 mL bottles of Body Spray, 1-75 mL bottle of Glaze Spray & 1-75 mL bottle of Glaze LT Spray.

Body Spray 75 mL

Specify:
Yellow (809-0252)
Yellow-Brown (809-0255)
Orange (809-0258)
Olive-Grey (809-0261)
Grey-Brown (809-0264)

Glaze Spray
(809-0265) 75 mL

Glaze Spray LT
(809-0713) 75 mL



VITA ENAMIC® Stains

VITA North America

Stains for reproducing natural shade nuances of hybrid ceramic restorations. Features a high level of individuality, which allows the anomalies and shade nuances of natural teeth to be individually characterized. Shade intensity can be controlled by the ratio of liquid to shade powder, ensuring the homogeneous application of shades.

(809-0301)Ea.

Contains: 6-2 Gm. stains (white, corn, brown, khaki, blue & black), 1-20 mL bottle of liquid, 1-20 mL bottle of glaze, 1 VITA VMLC cleaner, 20 disposable applicators, 2 brushes, 1 black porcelain mixing plate & 1 double-ended mixing spatula.

Stains Refill..... 2 Gm.

Specify:
White (809-0305) Khaki (809-0308)
Corn (809-0306) Blue (809-0309)
Brown (809-0307) Black (809-0310)

Accessories

Stains Liquid
(809-0302) 20 mL

VITA ENAMIC® Glaze

For surface sealing and increasing durability and brilliance.

(809-0300)Ea.
Contains: 1-20 mL bottle of glaze & 20 disposable applicators.

PROBLEM SOLVED.

Running a successful practice can feel overwhelming. You have a lot to worry about — providing the best patient care, increasing team productivity, and keeping up with advancements in dentistry — while also looking for ways to improve your business. Dentrrix has the answers. Our business tools solve real business problems and can help you to run a more profitable practice.

ONE DENTRIX

www.dentrrix.com



CEREC® Guide 3 Sleeve

Dentsply Sirona

Preassembled guide sleeve for the CEREC Guide 3 workflow. Enables guided insertion of an implant. After a CEREC Guide Bloc is milled, the sleeve is placed with an adhesive, eliminating the need for a drill key. Integrated depth stop ensures the precise transfer of an implant.

Sleeve..... Ea.

Specify:

Ankylos

ANK WD (811-8180)

ANK/XiVE ND (811-8178)

Astra Tech

AT EV ND (811-8181)

AT EV/XiVE WD (811-8182)

CEREC® Guide Drill Key Sets

Dentsply Sirona

Sets contain keys labeled with the manufacturer code that can only be used with the corresponding manufacturer and some universal drill keys. Universal drill keys are not labeled with a manufacturer code and are used in conjunction with the labeled keys in the key set. Made of grade 1.4301 stainless steel. Supplied unsterile.

Astra Tech: Facilitate (AT)

(587-2755) Ea.

Contains: 3 (S) drive keys, 1 of each diameter: 2mm, 3.2mm & 3.35mm; 6 (M) drive keys, 1 of each diameter: 2mm, 3.2mm, 3.35mm, 3.7mm, 3.85mm & 4.2mm; 8 (L) drive keys, 1 of each diameter: 2mm, 3.2mm, 3.35mm, 3.7mm, 3.85mm, 4.2mm, 4.7mm & 4.85mm; & 2 sterilizable boxes.

Biomet 3i: Navigator (B)

(587-2757) Ea.

Contains: 4 (S) drive keys, 1 of each diameter: 2mm, 2.75mm, 3mm & 3.25mm; 6 (L) drive keys, 1 of each diameter: 2mm, 2.75mm, 3mm, 3.25mm, 3.85mm & 4.25mm; 6 (M) drive keys, 1 of each diameter: 2mm, 2.75mm, 3mm, 3.25mm, 3.85mm & 4.25mm; & 2 sterilizable boxes.

Nobel Biocare (NB)

(587-2754) Ea.

Contains: 11 (L) drive keys, 1 of each diameter: 2mm, 2.8mm, 3mm, 3.2mm, 3.4mm, 3.6mm, 3.8mm, 4.1mm, 4.2mm, 4.6mm & 5mm; 5 (S) drive keys, 1 of each diameter: 2mm, 2.8mm, 3mm, 3.2mm & 3.4mm; 9 (M) drive keys, 1 of each diameter: 2mm, 2.8mm, 3mm, 3.2mm, 3.4mm, 3.6mm, 3.8mm, 4.1mm & 4.2mm; & 3 sterilizable boxes.

Straumann Guided Surgery Kit (ST)

(587-2752) Ea.

Contains: 3 (S) drive keys, 1 of each diameter: 2.2mm, 2.8mm, 3.5mm; 4 (L) drive keys, 1 of each diameter: 2.2mm, 2.8mm, 3.5mm, 4.2mm; 4 (M) drive keys, 1 of each diameter: 2.2mm, 2.8mm, 3.5mm, 4.2mm; & 2 sterilizable boxes.

Camlog Guide Drill Set (C)

Not suitable for use with CEREC Guide 1.

(587-2753) Ea.

Contains: 3 (L) drill keys in 4.5mm & 1 sterilizable box.

Scanbodies

Dentsply Sirona

Used with ScanPosts and TiBases for optical acquisition of implants. Enables the camera system to detect the position of the implant in order to design the anatomical features, size, and shape of a tooth restoration with computer software. Cap sits on top during imaging, and triangular tip dictates—to the software—the orientation of the implant. Available in connection sizes S and L; they correspond to the size of the hole, which is manufactured into the abutment blocks. Single-use only. Plastic.

White Caps

Compatible with the Bluecam & the following scanning systems: inEos X5, inEos Blue & CEREC AC with Bluecam.

White Caps36 per Pkg.

Specify:

Small (811-1036) Large (811-1035)

Grey Caps

Compatible with the CEREC AC Omnicam.

Grey Caps36 per Pkg.

Specify:

Small (811-1033) Large (811-1032)



Scanbodies

3shape

Implant scanning made easy. One thing is to take a digital implant impression instead of an analog one. Another is to use our intraoral scanner and scan bodies for it. It makes your digital implant workflow more efficient than ever. Use with original and 3rd party implant libraries. 100% titanium-made: autoclavable and reusable. Consistent performance and visible in x-rays.

Specify:

Biomet3i Certain Complete Kit (335-0333)

Biomet3i External Complete Kit (335-0334)

Camlog Complete Kit (335-0336)

Dentsply Implants (Friadent)

Astratech Kit (335-0337)

Dentsply Implants (Friadent)

Astratech EV Scanbody Kit (335-0338)

Dentsply Implants (Friadent)

XiVE Complete Kit (335-0340)

Dentsply Implants (Friadent)

Ankylos Scanbody Kit (335-0341)

Nobel Active Complete Kit (335-0343)

Nobel Replace Complete Kit (335-0344)

Osstem TS Complete Kit (335-0345)

Straumann Synocta Complete Kit (335-0346)

Zimmer TSV Complete Kit (335-0348)

ScanPosts

Dentsply Sirona

Impression post used intraorally to digitally capture the position of the implant in relation to the remaining teeth and soft tissue. Implant-specific scan post is screwed together with the implant, only for the purpose of optical detection.

Posts Ea.

Specify:

Ankylos

ANK

(587-4182)

Astra Tech EV

AT EV 3.0 S

(587-4185)

AT EV 3.6 S

(587-4186)

AT EV 4.2 L

(587-4187)

AT EV 4.8 L

(587-4189)

AT EV 5.4 L

(587-4191)

Astra Tech TX

AT OS 3.5/4.0 L

(587-2775)

AT OS 4.5/5.0 L

(587-2776)

BioHorizons

BH 3.0 S

(587-4177)

BH 3.5 L

(587-4178)

BH 4.5 L

(587-4179)

BH 5.7 L

(587-4181)

Biomet 3i

BC 3.4 S

(587-2788)

BC 4.1 L

(587-2789)

BC 5.0 L

(587-2790)

BO 3.4 L

(587-2777)

BO 4.1 L

(587-2778)

BO 5.0 L

(587-2780)

Dentsply Implants (Friadent)

FX 3.4 S

(587-2762)

FX 3.8 S

(587-2763)

FX 4.5 L

(587-2764)

FX 5.5 L

(587-2765)

Nobel Biocare

NB A 4.5 L

(587-2793)

NB A 5.0 L

(587-2794)

NB B 3.4 L

(587-2771)

NB B 4.1 L

(587-2774)

NBRS 3.5 L

(587-2766)

NBRS 4.3 L

(587-2767)

NBRS 5.0 L

(587-2769)

NBRS 6.0 L

(587-2770)

Osstem TS

O TS 3.5 L

(587-4175)

O TS 4.0 L

(587-4176)

Straumann

SBL 3.3 L

(587-2791)

SBL 4.1 L

(587-2792)

SSO 3.5 L

(587-2785)

SSO 4.8 L

(587-2786)

SSO 6.5 L

(587-2787)



FASTER ORDERING & MORE REWARDS

www.henryschein.ca

HENRY SCHEIN®
SOLUTIONS FOR HEALTH CARE PROFESSIONALS





TiBase Abutment System

Dentsply Sirona

Titanium base available for different implant systems and diameters. Includes a complete set with an abutment screw. Provides accurate digital acquisition of the implant position by the Scanbodies, extraorally on the model or intraorally, and abutment design using inLab or CEREC software.

Titanium Base..... Ea.

Specify:

Ankylos

ANK C GH1 S (587-4199)
ANK C GH2 S (587-4200)
ANK X GH1 S (587-4201)
ANK X GH2 S (587-4202)

Astra Tech EV

AT EV 3.0 GH1 S (587-4203)
AT EV 3.6 GH1 S (587-4204)
AT EV 4.2 GH1 L (587-4205)
AT EV 4.8 GH1 L (587-4206)
AT EV 5.4 GH1 L (587-4207)

Astra Tech TX

AT OS 3.5/4.0 L (587-2729)
AT OS 4.5/5.0 L (587-2730)

BioHorizons

BH 3.0 S (587-4195)
BH 3.5 L (587-4196)
BH 4.5 L (587-4197)
BH 5.7 L (587-4198)

Biomet 3i

BC 3.4 S (587-2743)
BC 4.1 L (587-2744)
BC 5.0 L (587-2745)
BO 3.4 L (587-2731)
BO 4.1 L (587-2732)
BO 5.0 L (587-2733)

Dentsply Implants (Friadent)

FX 3.4 S (587-2718)
FX 3.8 S (587-2719)
FX 4.5 L (587-2721)
FX 5.5 L (587-2722)

Medentika

MI 3.5/5.0 L (587-2749)

Nobel Biocare

NB A 4.5 L (587-2747)
NB A 5.0 L (587-2748)
NB B 3.4 L (587-2727)
NB B 4.1 L (587-2728)
NB RS 3.5 L (587-2723)
NB RS 4.3 L (587-2724)
NB RS 5.0 L (587-2725)
NB RS 6.0 L (587-2726)

Osstem TS

O TS 3.5 L (587-4192)
O TS 4.0 L (587-4194)

Straumann

S BL 3.3 L (587-2746)
S BL 4.1 L (587-2750)
SSO 3.5 L (587-2737)
SSO 6.5 L (587-2740)

Zimmer

ZTSV 3.5 L (587-2734)
ZTSV 4.5 L (587-2735)
ZTSV 5.7 L (587-2736)

FURNACES



CEREC SpeedFire

Dentsply Sirona

Compact, user-friendly furnace that can sinter and glaze full-zirconia restorations chairside. Induction technology allows short sintering and firing times. Can typically sinter a crown in 10–15 minutes and glaze in 9 minutes.

Complete Unit

(587-8299) Ea.



Programat® CS2

Ivoclar Vivadent

A glaze and crystallization furnace developed specifically for the dental office. Includes 20 individually programmable programs and 7 predefined programs, including the new speed program for the crystallization of IPS e.max® CAD in 15 minutes. Other features include a color touch display and a membrane-sealed keypad. 3-year warranty.

Complete Unit

(945-0509) Ea.

Contains: 1 furnace, 2 firing trays, 1 vacuum pump & 1 pair of tongs.

Automatic Temperature Checking Set

(945-8033) Ea.

Vacuum Pump

(945-0061) Ea.



Programat® CS4

Ivoclar Vivadent

Multipurpose furnace for the dental practice. Allows for the crystallization of IPS e.max® CAD and speed sintering of full-contour zirconia restorations, including IPS e.max® ZirCAD, in a single visit. Easy to use with a clearly arranged, membrane-sealed keypad. Graphic display ensures easy operation. Has a 30-program capacity: 8 preinstalled programs optimized for Ivoclar Vivadent materials and 22 open programs for individual programming. 3-year warranty.

Complete Unit

(945-1090) Ea.

Contains: 1 of each: CS4 furnace, power cord, temperature-checking set starter kit, IPS e.max speed-crystallization tray, sinter tray, sinter tray fork, "PrograBase 3" CD-ROM, USB cable, warranty card & operating instructions.

Sinter Tray

(945-1696) Ea.



VITA Smart.Fire®

VITA North America

Compact firing unit for chairside restorations. Indicated for crystallization, glaze, stain, and corrective firing of all common chairside materials. Preinstalled, material-specific programs for simple firing and 1-step touch-control operation. 2-year warranty.

Complete Unit

(809-0708) Ea.

Contains: 1 of each: firing unit, control unit & accessory set (2 platinum pins, 1 fibrous pad firing support, 1 firing tray, 1 firing plier & 1 cooling table).

CONTACT YOUR HENRY SCHEIN®
REPRESENTATIVE FOR ALL YOUR

CAD/CAM MATERIAL NEEDS

HENRY SCHEIN®
SOLUTIONS POUR LES PROFESSIONNELS DE LA SANTÉ

Rely on Us®

1-800-668-5558 • www.henryschein.ca





CEREC® SpeedPaste

Dentsply Sirona

Ceramic firing paste used to fix restorations on the restoration holder. Easy removal after firing and may be used with many ceramic materials. Not intended to be used for sintering.

Syringe

(587-2376)Ea.



Glaze Support

Dentsply Sirona

Supports restorations during the glazing process.

Single-Unit Restorations

(587-2382) Ea.

Multi-Unit Restorations

(587-2380) Ea.



GC Initial® Firing Foam

GC America

Refractory support material for ceramic firing.

Creamy consistency for easy application.

Dimensional stability during firing process. Easy to remove.

10 Gm. Syringe

(587-3701)Ea.



IPS e.max® Crystallization Trays and Pins

Ivoclar Vivadent

Silicon-nitride firing tray and pins for crystallizing IPS e.max® CAD restorations. Standard tray is for use on all Programat® furnace programs except the speed crystallization program. Small tray is used specifically with the speed crystallization program on all Programat® CS furnaces.

Standard Tray

(945-1229)Ea.

Contains: 1 tray & 9 pins.

Small Tray

(945-0585)Ea.

Pins

(945-9564) 18 per Pkg.



IPS Object Fix

Ivoclar Vivadent

Convenient substitute for firing pins or for use with veneer cases. Allows the restoration to be supported on the firing tray in a foam gel that promotes proper air flow around the restoration during the firing process.

IPS Object Fix

Convenient substitute for firing pins or for use with veneer cases. Allows the restoration to be supported on the firing tray in a foam gel that promotes proper air flow around the restoration during the firing process.

IPS Object Fix..... 10 mL

Specify:

Putty (945-1668)

Flow (945-8518)



Quick-Peg

VITA North America

High-heat insulator used to create custom firing pegs. Holds the restoration in place during firing and eliminates the use of hemostats. Indicated for use with all ceramic or PFM restorations, degassing, staining, and glazing.

10 mL Syringes

(809-0009)3 per Pkg.



Super Peg II

Harvest Dental

High-heat-resistant refractory material used for custom ceramic-firing stabilizing. Can be easily shaped and used with any type of porcelain. Is stable up to 1300°C (2375°F), which eliminates firing distortion and contamination.

12cc Syringes

(671-0079) 4 per Pkg.



STAY CONNECTED with

HENRY SCHEIN®

Stay informed on the most up-to-date news in the dental industry with our many social network platforms. Just one more reason to **Rely on Henry Schein.**

[www.Facebook.com/HenryScheinCanada](https://www.facebook.com/HenryScheinCanada)

[@hscinecanada](https://www.instagram.com/hscinecanada)

[Henry Schein Canada](https://www.youtube.com/HenryScheinCanada)

www.linkedin.com/company/HenryScheinCanada





Chairside Dental Designer Set

Kerr Rotary

Offers dental professional with an effective and efficient method of finishing high-quality CAD/CAM restorations (either composite or ceramic) chairside. Kit is delivered in a sterilizable, aluminum bur block.

LS-7576

(228-0927)Ea.



Two Striper® Universal All-Ceramic Prep Kit

Premier Dental

Assorted kit that contains all of the instruments needed to systematically create fast and accurate preparations for chairside CAD/CAM or lab-fabricated ceramics.

Complete Kit

(378-3325)Ea.

Contains: 1 of each diamond: 2014.5C, 703.8KRC, 703.8KRF, 770.8C, 770.8F, 263.8F, 200.3F, 285.5F & 292.3F.



Universal CAD/CAM Prep & Finishing Kit

SS White

Adds efficiency and organization to CAD/CAM preparations with a streamlined selection of rotary instruments for milled ceramic crowns, inlays, and partial coverage restorations. Kit features the innovative SS White® Great White® 1.5mm- and 2mm-depth cutter, along with selected SS White® carbide burs, diamond instruments, and Jazz® polishers, organized in an autoclavable metal bur block.

#64502 Complete Kit

(999-1765)Ea.

Contains: 1 of each bur or polisher: #5040 GWDC1, #15041 GWDC2, #16884 FG #169, #15675-5 finishing bur #7675, #18193-5 Great White® Ultra #379-023, #18182-5 Great White® Ultra #856-016S, #881-016VF diamond, #89012 Jazz® polisher (fine flame P3S), #89008 Jazz® polisher (medium flame P3S) & #89004 Jazz® polisher (coarse flame P3S).



Maestro® Ceramic Abrasives Kit

HENRY SCHEIN®

High-performance, synthetically bound diamond abrasives for removing metal. Specially designed for zirconia and alumina restorations without generating heat. Coarse abrasives for gross reduction and medium abrasives for finishing.

Complete Kit

(112-5793)Ea.

Contains: 1 of each: coarse disc, coarse cone, medium disc & medium cone.

Disc Refills..... 3 per Pkg.

Specify:

Coarse (112-5787) Fine (112-5790)

Cone Refills..... 3 per Pkg.

Specify:

Coarse (112-5789) Fine (112-5791)



The ceredoctors.com Prep Kit

Meisinger

Developed in conjunction with ceredoctors.com, this kit is intended for use with lithium-disilicate, silicate, leucite-reinforced, and feldspathic ceramics. All instruments are labeled and organized for easy reference in an anodized bur block.

SC02 Prep Kit

Contains diamonds for inlays, onlays, and full-coverage crowns in a convenient arrangement for short processing times.

(701-3902)Ea.

Contains: 1 of each FG bur (856F-016, 856F-018, 856F-021, 856H-016, 856H-018, 801F-018, 801F-023, 845RH-018, 845RH-025, 845RF-018, 828Y-017, 833-023, 833F-023, 828G-017, 850F-012, 850H-012 & 850H-021).



Occlusal Reduction Kit

Meisinger

Kit assortment developed by Athas N. Kometas, DMD. Allows for precise occlusal reductions. Occlusal reduction burs are available in 1.5mm, 1.8mm, 2.0mm, 2.2mm, and 2.4mm depths that correspond to proper occlusal reduction for various restorative materials. A 1.0mm bur is also included for labial reduction.

2565 Reduction Kit

(701-0546)Ea.

Contains: 1 of each: 828G-017-FG (1.5mm), 828O-017-FG (1.8mm), 828Y-017-FG (2.0mm), 828B-017-FG (2.2mm), 828R-017-FG (2.4mm) & 828W-017-FG (1.0mm).



The Klim Institute Prep and Finishing Kit

Meisinger

Kit assortment developed by James Klim, DDS. Provides efficient shaping with Black Cobra diamonds and corresponding complimentary-shaped carbide finishing burs for lathe-like finished margins and line angles.

JK02 Complete Kit

(701-1573)Ea.

Contains: 1 of each: 833F-023-FG, 833C-023-FG, 898F-014-FG, 898C-014-FG, HM44E-023-FG, HM246-009-FG, 849G-016-FG, HM375R-018-FG, HM375R-016-FG, HM48L-012-FG, 909G-040-FG, 828TY-017-FG, HMB21R-012-FG, B855-018-FG, B850-016-FG, B862-014-FG, B830-023-FG & 828TG-017-FG.

**It's Quick,
It's Easy,
It's Accurate**



ZAHN CANADA®
A DIVISION OF HENRY SCHEIN®

1-800-496-9500 • www.zahncanada.ca

PROBLEM SOLVED.

Running a successful practice can feel overwhelming. You have a lot to worry about — providing the best patient care, increasing team productivity, and keeping up with advancements in dentistry — while also looking for ways to improve your business. Dentrrix has the answers. Our business tools solve real business problems and can help you to run a more profitable practice.

 **ONE DENTRIX**
www.dentrrix.com

Alpen® ShapeGuard™ Zirconia Plus System

Coltene

Two-step zirconia polishing system for fast, easy, and intuitive finishing and polishing. Flexible polishing lamellas adapt to any surface for optimal polishing. RA polishers come in a 1/10mm diameter in both medium and fine grits. Able to be sterilized and reused, aiding in a higher service life span.

Complete System
(129-4344)Ea.

Contains: 4-#140 polishers (2 of each: Step 1-Medium & Step 2-Fine).

Zirconia Refills 4 per Pkg.

Specify:
Step 1-Medium (129-4342)
Step 2-Fine (129-4343)

Diamond Silicone Chairside Polishing Kit

Dedeco

Specially designed for removing scratches from porcelain, zirconium, and e.max® and polishing to a high-glaze finish.

1183 Complete Kit
(142-1042)Ea.

Contains: 9 polishers (1 of each: 1107, 1108, 1109, 1127, 1128, 1129, 1147, 1148 & 1149) & 1-12-hole Lucite bur block.



Diamond Silicone Extraoral Polishing Kit

Dedeco

A 3-step, diamond-filled, silicone rubber system for adjusting, smoothing, and high-glaze polishing of porcelain, zirconium, lithium disilicate, and e.max®. Eliminates the need for polishing paste or reglazing.

1181 Complete Kit

(142-0186)Ea.

Contains: 6 polishers (1 of each: 1103, 1105, 1123, 1125, 1143 & 1145).



eZr™ Zirconia Polishing Kit

Garrison Dental Solutions

Three-step diamond polishing system for processing zirconium oxide and lithium disilicate. Grinders and polishers are ideal for adjusting and correcting chairside restorations such as inlays, onlays, veneers, and crowns. eZr™ has a high material-removal rate and produces a high-gloss finish on zirconia without polishing paste.

Complete Kit

(300-0264)Ea.

Contains: 1 of each: coarse disc grinder, coarse flat end taper grinder, medium large point, medium disc, medium cup, fine large point, fine disc & fine cup).

Refills

Blue Coarse Grinders

Disc

(300-0312) 5 per Pkg.

Flat-End Taper

(300-0290) 5 per Pkg.

Assorted Shapes..... 5 per Pkg.

Specify:

Large Point (300-0270) Cup (300-0291)

Disc (300-0307)

Assorted Shapes..... 5 per Pkg.

Specify:

Large Point (300-0289) Cup (300-0292)



OptraFine® Diamond Polishing System

Ivoclar Vivadent

High-performance diamond-polishing system. Designed specifically for ceramic materials. Creates a smooth surface and provides an aesthetic high-gloss finish. Autoclavable and reusable.

Assortment Pack

(945-0035)Ea.

Contains: 2 of each finisher: cup, flame & disc; 2 of each polisher; cup, flame & disc; 6 HP brushes & 1-1 mL polishing paste.

Refills..... 10 per Box

Specify:

Finishers

Cup (585-9755)

Disc (585-9756)

Flame (585-9757)

Polishers

Cup (585-9758)

Disc (585-9759)

Flame (585-9760)

HP Brush Nylon Refill

(585-9761) 10 per Box

High Polishing Paste Refill

(585-9762) 1 mL



Zir-Cut™ Polishers

Kerr Rotary

A 2-step system that removes scratches, smooths rough surfaces, and creates a high-gloss surface on all zirconia crowns, bridges, substructures, and abutments. Unique diamond-abrasive mixture brings zirconia materials to a high shine.

LS-517 Zirconia Polisher Set

(228-1975)Ea.

Contains: 3 polishers (1 of each: flame, cup & disc).

Refills..... 3 per Pkg.

Specify:

Prepolishing – Green

Cup (228-1980)

Disc (228-1982)

Flame (228-1978)

High Shine – Orange

Cup (228-1989)

Flame (228-1984)

Disc (228-1991)

The CAD/CAM Finishing Kit

Meisinger

This universal contouring and finishing kit was created to provide a wide range of solutions for today's ceramic restorations. The kit provides all of the necessary tools to produce excellent results for both anterior and posterior crowns.

(701-3906) Ea.

Contains: 1 of each: 9770G-260-HP, 9770M-260-HP, 9770F-260-HP, 9736G-150-HP, 9613V-220-UNM, 303RF-050-HP, 9734G-040-HP, 882L-018-HP, 850-018-HP, 943DC-100-HP, 9771G-170-HP, 9771M-170-HP, 9771F-170-HP, 9771C-170-HP, 863-016-HP, 805S-018-HP & 936SC-200-HP.



LUSTER® for Silicate Ceramics

Meisinger

Polishing system developed for use on post-crystallized lithium-disilicate restorations. Diamond polishers are suitable for intraoral and extraoral use.

LUS100 Intraoral Adjusting & Polishing Set
(701-3448) Ea.

Contains: 1 of each: DP002, DP007, DP005, DP008, Z0801, Z850F & Z833F.

LUS80 Extraoral HP Laboratory Set
(701-3340)Ea.

Contains: 1 of each: 9736G HP, DPO02 HP, DPO03 HP, DPO05 HP & DPO06 HP.



The Klim Institute CAD/CAM Lab Kit

Meisinger

Kit developed for finishing CEREC® restorations, according to James Klim, DDS. Intended for use with lithium-disilicate and silicate ceramics, Leucite®-reinforced ceramics, and feldspathic ceramics.

(701-3459) Ea. **509.99**

Contains: 1 of each: 110-190-HP, 805S-018-HP, 850-018-HP, 882L-018-HP, 936SC-200-HP, 943DC-100-HP, 9734G-040-HP, 9735G-050-HP, 9736G-150-HP, 9770G-260-HP, 9770M-260-HP, 9771G-170-HP, 9771M-170-HP, 9771F-170-HP & 9771C-170-HP.

LUSTER® for Zirconia

Meisinger

Developed for use on layered and monolithic zirconium-oxide restorations. Performs tasks quickly without creating excessive heat. Two-stage polishers can be used without polishing paste. Intended for use with zirconium.

#LUS91, Intraoral Complete Kit
(701-3447) Ea.

Contains: 1 of each: Z0801 FG, Z0850 FG, Z0833 FG, DCA04 RA, DCA13 RA, DCA10 RA & DCA14 RA.

LUS85, Extraoral HP Laboratory Set
(701-3355) Ea.

Contains: 1 of each: 9736H HP, Z652R HP, DCA04 HP, DCA06 HP, DCA10 HP & DCA12 HP.

LUSTER® The Klim Institute

Meisinger

Kit developed according to James Klim, DDS. Available in 3 shapes and 4 grits for precise adjusting and polishing. Creates a long-lasting final shine. Intended for use with porcelain and ceramic-like materials.

JK01 Complete Kit
(701-3467) Ea.

Contains: 1 of each: 9771G-140-RA, 9744G-060-RA-GR/O, 9746G-100-RA-GR/O, 9771F-140-RA, 9744F-060-RA-R/O, 9746F-100-RA-R/O, 9771M-140-RA, 9744M-060-RA-BL/O, 9746M-100-RABL/O, 9771C-140-RA, 9744C-060-RA-Y/O & 9746C-100-RA.

LUSTER® Intraoral Polishing Kit

Meisinger

Diamond-impregnated 4-step polishing system. For polishing high-gloss ceramic restorations in regions difficult to access. Flexible disks contour and adjust to every dip and curve of the surface but still preserve its original structure.

LUS06 Intraoral Kit
(701-3862) Ea.

Contains: 4 RA polishers (1 of each: Pre-polisher, Polisher, Final High-Shine Polisher & Mirror-Finish Polisher) & 1 bur holder.

Intraoral Refills.....2 per Pkg.

Specify:
Extra Fine-Yellow (701-3904)
Fine-Red (701-3911)
Medium-Blue (701-3908)
Coarse-Green (701-3903)

LUS05 Extraoral Kit
(701-3861)Ea.

Contains: 4 RA polishers (1 of each: Pre-polisher, Polisher, Final High-Shine Polisher & Mirror-Finish Polisher) & 1 bur holder.

Extraoral Refills.....2 per Pkg.

Specify:
Extra Fine-Yellow (701-3850)
Fine-Red (701-3859)
Medium-Blue (701-3848)
Coarse-Green (701-3860)



Twist Polishers Kits for Hybrid Ceramics

Meisinger

Developed according to Greg Campbell, DDS. Provides efficient processing of poly ceramics and brings out the luster of restorations. Intended for use with poly ceramics and flowable composites.

GC05 Extraoral Kit
(701-3914) Ea.

Contains: 1 of each: 9735G-050-HP-P, 9769F-170-HP, 9787-145-HP, 805S-018-HP & 9769M-170-HP.

GC06 Intraoral Kit
(701-3915)Ea.

Contains: 1 of each: 9762M-040-RA, 9769M-140-RA, 9769F-140-RA, 9762F-040-RA-GREY, 805F-014-FG, 862F-010-FG & 830F-021-FG.



Twist Polishers Kits for Silicate Ceramics

Meisinger

Developed in collaboration with Greg Campbell, DDS. Provides easy and efficient processing of silicate ceramics and can be used on precrystallized and crystallized restorations. Intended for use with silicate ceramics, Leucite®-reinforced ceramics, and feldspathic ceramics.

GC01 Extraoral Kit
(701-3909)Ea.

Contains: 1 of each: 9735G-050-HP-P, DPO10-170-HP, 9771C-170-HP, 805S-018-HP, DPO09-170-HP & 9771M-170-HP.

GC02 Intraoral Kit
(701-3910) Ea.

Contains: 1 of each: DPO02-040-RA-R/O, 9771M-140-RA, 9771C-140-RA, DPO05-040-RA-G/O, 805F-014-FG, 862F-010-FG & 830F-021-FG.



Twist Polisher Kits for Zirconia

Meisinger

Developed according to Greg Campbell, DDS. Provides efficient processing of zirconia and high-gloss results. Intended for use with monolithic zirconium ceramics and zirconium-reinforced lithium-silicate (ZLS) materials.

GC03 Extraoral Kit
(701-3912)Ea.

Contains: 1 of each: 9735H-050-HP, DCA12-170-HP, 9771C-170-HP, 805S-018-HP, DCA06-170-HP & 9771M-170-HP.

GC04 Intraoral Kit
(701-3913)Ea.

Contains: 1 of each: DCA04-040-RA-BL/O, 9771M-140-RA, 9771C-140-RA, DCA10-040-RA-BRN/O, 805F-014-FG, 862F-010-FG & 830F-021-FG.



Poli-Pro Disks™ Composite Polishing System



Premier Dental

Safe, nonmetal center disks that will not gouge or discolor the restoration. Wider abrasive surface allows for faster polishing. Thin and flexible disks conform to tooth anatomy. Single-use design avoids cross contamination.

System Kit

(378-6953)Ea.

Contains: 180 disks (20 of each size) & 4 RA mandrels.

Refills.....40 per Box

Specify:

Coarse-Black

1/2" Grit Up (378-1107)

5/16" Grit Up (378-6395)

Medium-Blue

1/2" Grit Up (378-4930)

1/2" Grit Down (378-1915)

5/16" Grit Up (378-0785)

Fine-Red

1/2" Grit Double-Sided (378-4424)

5/16" Grit Double-Sided (378-7761)

X-fine-White

1/2" Grit Double-Sided (378-5014)

5/16" Grit Double-Sided (378-0331)

RA Mandrels

(378-2184) 4 per Pkg.

FG Mandrels

(378-1202) 4 per Pkg.

JAZZ™ P3S Porcelain, Metal 3-Step Polishing System



SS White

Diamond impregnated to achieve brilliant results on porcelain, metal, and restorations with minimal pressure. Eliminates the need for multiple polishing systems. Reusable for multipatient use.

#89016 Assortment Kit

(999-0671)Ea.

Contains: 12 polishers (1 of each shape & grit).

Refills.....3 per Pkg.

Specify:

Coarse-Green

Cup (999-0674)

Knife (999-0675)

Medium-Brown

Cup (999-0678)

Knife (999-0679)

Fine-White

Flame (999-0681)

Cup (999-0682)

Knife (999-0683)

JAZZ® Polishers ZA2S System



SS White

Specially formulated polishers for full-contour zirconia restorations, featuring a high diamond grit for more efficiency and durability. Achieve a high gloss and eliminate the need for glazing.

Sterile..... 3 per Pkg.

Specify:

Medium

Flame FG (999-2590)

Small Flame-RA (999-2594)

Fine

Flame-FG (999-2591)

Small Flame-RA (999-2595)

Sterile..... 3 per Pkg.

Specify:

Medium

Flame-RA (999-2627)

Cup-RA (999-2626)

Knife Edge-RA (999-2622)

Fine

Flame-RA (999-2628)

Cup-RA (999-2623)

Knife Edge-RA (999-2625)

Non-Sterile3 per Pkg.

Specify:

Medium

Flame-FG (999-2613)

Small Flame-RA (999-2614)

Fine

Flame-FG (999-2612)

Small Flame-RA (999-2611)

Non-Sterile3 per Pkg.

Specify:

Medium

Flame-RA (999-2619)

Cup-RA (999-2617)

Knife Edge-RA (999-2615)



JAZZ® Pro Brush Polishers PMC2S

SS White

Diamond-impregnated brushes for polishing CAD/CAM, porcelain, composites, and metals. Produce glossy surfaces after minimal processing time. Low heat generation with minimal polishing paste required. RA latch-type.

Sterile Medium 3 per Pkg.

Specify:

Small Cup (999-2598) Flame (999-2597)

Large Cup (999-2596)

Sterile Fine..... 3 per Pkg.

Specify:

Small Cup (999-2587) Flame (999-2589)

Large Cup (999-2588)

Diamond Twist SCL™



Premier Dental

Designed for extraoral prepoly smoothing and glaze-like final polish on castable ceramics, CAD/CAM porcelain restorations, or composite crowns requiring multisurface polishing.

Complete Kit

(378-1266) Ea.

Contains: 1-6 Gm. twist-stick Diamond Twist polishing paste, 2 HP mandrels, 4 prepolishing fibra shapes & 4 felt paste applicators.

SCL Polishing Paste Refill

(378-4013)6 Gm.

Diamond Twist SCO™



Premier Dental

Designed for intraoral polishing of indirect porcelain restorations and direct nanohybrid composites.

Complete Kit

(378-7024)Ea.

Contains: 1-3 Gm. syringe Diamond Twist polishing paste, 1 RA mandrel, 6 polish applicator brushes & 6 Micro-Cloth disk paste applicators.

SCO Polishing Paste Refill

(378-8416)3 Gm.

POLISHING PASTE

DiaSheen™ Diamond Polishing Paste



Kerr Rotary

Formulated to create a high gloss on composites and porcelain. Firm consistency prevents splatter and waste. Ideal for use with stiff bristle or felt wheels.

3 Grits.....3 Gm. Jar

Specify:

Fine (228-9542)

Fine-Soft (228-4428)

Medium (228-4284)

HENRY SCHEIN®
BRAND PRODUCTS



HAVE YOU SWITCHED TO
**HENRY SCHEIN
BRAND**
PRODUCTS YET?

1-800-668-5558 | henryschein.ca

Rely on Us™

©2021 Ivoclar Vivadent, Inc. IPS e.max and Ivoclar Vivadent are registered trademarks of Ivoclar Vivadent, Inc. Rev. 2
* See Official Rules at https://www.ivoclarvivadent.com/NA/makeitemax_contest

Make it e.max.®

Dr. Kristine Aadland, DMD Vancouver, WA

"I love the options e.max provides in the different translucencies and shades. I can also sleep better at night knowing I have delivered a beautiful smile that will last. I have been using e.max for most of my dental career and never regretted it, nor have I found a truly universal product that I love more. I believe in e.max through and through."

Learn more @ [makeitemax.com](https://www.makeitemax.com)



You could WIN a Peloton Bike+

Plus other prizes!* Post a video or photo saying why you #makeitemax on FB, IG or Twitter. Use #makeitemax #ivoclarvivadentna

[ivoclarvivadent.com](https://www.ivoclarvivadent.com)
Making People Smile

ivoclar
vivadent®



Calibra® Ceram
Dentsply Sirona Restorative

Adhesive resin cement. Light cured, dual cured, or self-cured paste system that does not need a self-cure activator. Shade Stable™ Technology. Wide tack-cure window up to 10 seconds and extended 45-second gel phase for stress-free cleanup. Designed for bonding and luting to indirect restorative materials including ceramic, composite, and metal-based inlays, onlays, crowns, bridges, and posts.

Combination Kit

(222-0872)Ea.

Contains: 1-4.5 Gm. automix syringe, 1-5 mL bottle of Prime&Bond Elect® adhesive, 25 flocked applicator tips, 10 mixing tips & 1 dispensing well.

Refill PackageEa.

Specify:

Light (222-0873)

Medium (222-0874)

Translucent (222-0875)

Opaque (222-0878)

Bleach (222-0879)

Contains: 1-4.5 Gm. automix syringe & 10 mixing tips.



Calibra® Universal
Dentsply Sirona Restorative

Dual-cure, self-adhesive resin cement with fluoride. Shade Stable™ technology helps avoid shade shifting. Suitable for permanent cementation of metal, PFM, resin/composite, ceramic and porcelain inlays, onlays, crowns and bridges, and endodontic posts. Wide tack-cure window up to 10 seconds. Extended 45-second gel phase for excess cement cleanup.

Refill Package Ea.

Specify:

Light (222-0797)

Medium (222-0802)

Translucent (222-0803)

Opaque (222-0809)

Bleach (222-0848)

Contains: 2-4.5 Gm. automix syringes & 20 mixing tips.



G-CEM LinkAce™

GC America

Self-adhesive resin cement. High bond durability to zirconia in one step. Easy removal of excess cement. Can be used on all ceramic, metal, and indirect composite restorations. Dual cure.

Handheld Syringe RefillsEa.

Specify:

A2 (333-3390)

AO3 (333-3392)

BO1 (333-3393)

Translucent (333-3391)

Contains: 2-4.6 Gm. syringes, 15 regular automix mixing tips & 5 endo tips.



G-CEM LinkForce™

GC America

Dual-cure adhesive resin cement. Strong, aesthetic, and universal solution for all adhesive cementations. Dual-cure activator offers high bond strength through light- and self-cure modes. Superior wear resistance, easy cleanup, and thin film thickness. Tooth-like fluorescence. Indicated for all types of ceramic, resin, and metal-based inlays, onlays, and crown and bridge; composite veneers (up to 2 teeth); and permanent cementation of crowns and bridges on implant abutments.

System Kit

(333-4274)Ea.

Contains: 3-5 mL syringes (1 of each shade: A2, translucent & opaque), 45 automix tips, 5 automix endo tips mixed with intraoral tips, 1-5 mL bottle of G-Premio Bond, 1-5 mL bottle of G-Multi Primer, 1-3 mL G-Premio Bond DCA, 1-3.6 mL GC etchant syringe, 25 metal dispensing tips with 2 covers, 3-1.2 mL try-in paste syringes (1 of each shade: A2, translucent & opaque), 50 disposable applicators & 20 disposable dispensing dishes.

Starter KitEa.

Specify:

A2 (333-4269) Translucent (333-4271)

Contains: 1-5 mL syringe, 1-5 mL bottle of G-Premio Bond, 1-5 mL bottle of G-Multi Primer, 20 automix tips, 50 disposable applicators & 20 disposable dispensing dishes.

Refill PackageEa.

Specify:

A2 (333-4253)

Translucent (333-4254)

Opaque (333-4255)

Bleach (333-4256)

Contains: 1-5 mL syringe & 20 automix tips.

G-Multi Primer-Bottle (333-4267) 5 mL

G-Premio BOND Bottle Refill (333-4938)5 mL



Maxcem™ Elite™ Chroma

Kerr Restoratives

Dual-cure, self-etch, self-adhesive resin cement. Indicator tells the optimal time to clean up: when the pink color fades at the gel state. One-peel™ feature ensures easier removal of excess cement. Compatible with all indirect restorations including ceramics, PFMs, metal restorations, and CAD/CAM materials. Radiopaque.

Syringe Refills Ea.

Specify:

Clear (124-5437)

Yellow (124-5440)

White (124-5439)

Contains: 2-5 Gm. syringes, 16 pink automix tips (8 regular & 8 wide) & 8 intraoral tips.

Bulk Kit

(124-5443)Ea.

Contains: 4-5 Gm. syringes, 24 pink automix tips (16 regular & 8 wide), 8 intraoral tips & 8 root canal tips.



Meron Plus QM

VOCO

Resin-modified glass-ionomer cement. Self-adhesive material with enhanced bond strength to dentin. Moisture and acid resistant. Bonds without primer or conditioner. Indicated for cementing metal-based inlays, onlays, high-strength (zirconium dioxide-based) all-ceramic crowns and bridges, all-ceramic inlays, metal, ceramic, and fiber posts. Automix delivery system.

Quickmix Syringe

(999-2543)Ea.

Contains: 1-8.5 Gm. syringe & 10 type-10 mixing tips.

(999-2544) 3 per Pkg.

Contains: 3-8.5 Gm. syringes & 30 type-10 mixing tips.

**Multilink®
Automix Next
Generation
Universal Resin Cement**

Ivoclar Vivadent

Dual-cure, self-etch, universal resin cement. Provides the security of high bond strength with a fast, simple application technique. Automix delivery provides accurate application with no mixing. Indicated for cementing metal, metal-free, and zirconia crowns as well as inlays, onlays, bridges, and posts.

NG System PackageEa.

Specify:
Transparent (945-0462)
Yellow (945-0463)
White (945-0460)
Opaque (945-0464)

Contains: 1-9 Gm. syringe, 1-3 Gm. bottle of each primer (A & B), 1-5 mL bottle of Monobond Plus, 1-2.5 Gm. liquid strip, 5 liquid strip application cannulas, 15 tapered mixing tips, 5 root-canal mixing tips, 50 microbrushes (25 small & 25 large), 1 mixing plate, 2 Opra Sticks & 1 mixing pad.

NG Refill Syringe.....Ea.

Specify:
Transparent (945-0466)
Yellow (945-0467)
White (945-0465)
Opaque (945-0468)

Contains: 1-9 Gm. syringe & 15 automix mixing tips.

Monobond Plus

(945-0200)5 Gm.



**NX3® Nexus®
Third Generation**

Kerr Restoratives

Universal, adhesive resin cement. High bond strengths and wear resistance. For all indirect applications, including veneers. Can be used with a total-etch or self-etch bonding agent.

5 Gm. Dual-Cure Syringe Ea.

Specify:
Clear (123-5340)
Yellow (123-8385)
White (123-3907)
Bleach (123-0962)
White Opaque (123-0964)



ZAHN CANADA®
A DIVISION OF HENRY SCHEIN®

1-800-496-9500 • www.zahncanada.ca

Panavia™ V5

Kuraray America

Dual-cure, universal adhesive resin cement. Fluoride releasing and easy cleanup. MDP-based Panavia V5 Tooth Primer provides high bond strength for all tooth structures in self-cure mode. Color-stable shades.

Standard KitEa.

Specify:
Universal (A2) (721-0667)
Clear (721-0668)

Contains: 1-4.6 mL syringe of Panavia V5 Paste, 1-1.8 mL syringe of Panavia V5 Try-In Paste, 1-4.0 mL bottle of Panavia V5 Tooth Primer, 1-4.0 mL bottle of Clearfil Ceramic Primer Plus, 1-3.0 mL syringe of K-Etchant, 15 mixing tips, 5 endo tips, 50 applicator brushes, 20 needle tips & 1 mixing dish.

Introductory KitEa.

Specify:
Universal (A2) (721-0669)
Clear (721-0670)

Contains: 1-2.4 mL syringe of Panavia V5 Paste, 1-2.0 mL bottle of Panavia V5 Tooth Primer, 1-2.0 mL bottle of Clearfil Ceramic Primer Plus, 10 mixing tips, 50 applicator brushes & 1 mixing dish.

Refill Package Ea.

Specify:
Universal (A2) (721-0656)
Clear (721-0657)
White (721-0659)
Opaque (721-0661)

Contains: 1-4.6 mL syringe of Panavia V5 Paste & 20 mixing tips.

Try-In Paste

1.8 mL Syringe Ea.

Specify:
Universal (A2) (721-0662)
Clear (721-0663)
White (721-0665)

Contains: 1-4.6 mL syringe & 20 mixing tips.

Clearfil™ Ceramic Primer Plus-Bottle
(721-0653)4 mL

Panavia™ V5 Tooth Primer-Bottle
(721-0654)4.0 mL



3M

**RelyX™
Luting Plus Automix**

Resin-modified, glass-ionomer cement packaged in a handheld automix syringe. Tack light-cure option for faster removal of excess cement.

Trial Kit #3535TK

(777-0455) Ea.

Contains: 1-8.5 Gm. handheld automix syringe & 12 automix mix tips.

Value Package #3535

(777-0454) Ea.

Contains: 3-8.5 Gm. handheld automix syringes & 36 mixing tips.

3M

RelyX™ Unicem Maxicap™

Self-adhesive, universal resin cement in larger capsules, suited for multiple-unit restorations, bridgework, and single-unit restorations that require a longer working time.

Introductory Kit #56830

(378-8407) Ea.

Contains: 10-936 mg. capsules of A2, 6-936 mg. capsules of translucent, 4-936-mg. capsules of A3, 1 activator & 1 applier.

Maxicap™ Refill Ea.

Specify:
A2 Universal #56834 (378-5650)
A3 Opaque #56833 (378-2232)
Translucent #56832 (378-5663)
Assorted #56831 (378-4230)

Contains: 20-936 mg. capsules.



3M

RelyX™ Unicem Aplicap™

Dual-cure, self-adhesive, universal resin cement. Can be used for virtually every indication: metal, ceramic, composite material, crowns & bridges, inlays/onlays, posts, and screws. Easy to handle. No pretreatment or absolute moisture control is necessary. Combines the advantages of both conventional cements (ease of use) and resin cements (performance).

Aplicap Trial Kit #56829

(777-2940) Ea.

Contains: 10-295 mg. (A2), 6-295 mg. (Translucent), 4-295 mg. (A3e), 1 activator, 1 applier instrument set, 1 technique guide & instructions.

Aplicap Introductory Kit #56814

(777-8478)Ea.

Contains: 24-295 mg. (A2), 16-295 mg. (Translucent), 10-295 mg. (A3), 1 activator, 1 applier instrument set, 1 technique card & instructions.

Aplicap™ Refills50 per Box

Specify:
A2 Universal #56818 (378-3757)
A3 Opaque #56817 (378-3143)
Translucent #56816 (378-8676)
Assorted #56815 (378-8736)

Contains: 50-295 mg. capsules.

Aplicap™ Refills20 per Box

Specify:
A1 #56819 (378-9981)
White Opaque #56820 (378-1329)

Contains: 20-295 mg. capsules.

ORDER ONLINE!



www.henryschein.ca HENRY SCHEIN® SOLUTIONS FOR HEALTH CARE PROFESSIONALS



3M RelyX™ Luting Plus

Resin-modified, glass-ionomer cement. Packaged in the convenient Clicker™ dispenser. Easy to dispense, mix, and load. Chemically bonds to tooth and restoration. Tack light-cure option for faster removal of excess cement.

Trial Kit #3525TK

(777-2737)Ea.

Contains: 1–11.5 Gm. Clicker dispenser of cement, technique guide & set of instructions.

Refill Package #3525

(777-6701)Ea.

Contains: 2–11.5 Gm. Clicker dispensers of cement & set of instructions.

3M RelyX™ Unicem 2 Automix

Self-adhesive resin cement. One-step automix material in an easy-to-use handheld syringe. Delivers high adhesion, resistance to discoloration, and long-term stability. Offers enhanced bond strength and aesthetics and easy removal of excess cement.

Starter Kit #56850

(378-2957)Ea.

Contains: 1–5 mL syringe of shade A2, 5 regular mixing tips, 10 wide mixing tips, 10 endo tips, 5 intraoral tips, 20 fiber post refills (5 of each size: 0, 1, 2 & 3) & 5 drills (1 of each size: 0, 1, 2, 3 & universal).

Refill.....Ea.

Specify:

A2 #56846 (378-2966)

A3 #56847 (378-2981)

Translucent #56848 (378-2962)

Contains: 1–5 mL syringe, 10 regular mixing tips, 5 wide mixing tips & 5 endo tips.

Value Pack.....Ea.

Specify:

A2 #56849 (378-2951)

A3 #56857 (378-2960)

Translucent #56858 (378-2946)

Contains: 3–5 mL syringes, 30 regular mixing tips, 20 wide mixing tips, 15 endo tips & 5 intraoral tips.



New!

RelyX™ Universal

Universal resin cement.

Works both as a standalone, self-adhesive cement and as an adhesive cement when combined with 3M™ Scotchbond™ Universal Plus Adhesive. Together, they eliminate the hassle and confusion of multiple resin cements, primers, and adhesives to simplify direct and indirect restorations. Excellent bond strength, high esthetics, and easy excess cleanup.

Trial Kit.....Ea.

Specify:

Translucent #56969 (378-0958)

A1 #56970 (378-0959)

Contains: 1–3.4 Gm. syringe, 1–3 mL Scotchbond™

Universal etch gel, 1–1.5 mL Scotchbond Universal

adhesive, 25 blue etchant tips, 15 Universal Micro

Mixing tips, 5 elongation tips & 50 orange microbrushes.

Refill.....Ea.

Specify:

Translucent #56971 (378-0951)

White Opaque #56974 (378-0954)

A1 #56972 (378-0950)

A3/Opaque #56973 (378-0953)

Contains: 1–3.4 Gm. syringe, 15 Universal Micro

Mixing tips & 5 elongation tips.

Value Pack.....Ea.

Specify:

Translucent #56977 (378-0955)

A1 #56978 (378-0957)

Contains: 3–3.4 Gm. syringes, 45 Universal Micro Mixing tips & 15 elongation tips.



3M

RelyX™ Unicem 2 in Clicker™ Dispenser

Self-adhesive universal resin cement in the convenient Clicker™ dispenser. No need to pretreat etching, priming, or bonding. Provides a convenient delivery system and a reliable 1:1 mix ratio.

Complete Package

11 Gm.Ea.

Specify:

A2 #56875 (378-0300)

A3 #56876 (378-0301)

Translucent #56874 (378-0299)

3M RelyX™ Veneer Cement

Permanent, light-cure-only resin luting material for cementing ceramic or composite veneers. Formulated for accurate color matching of try-in paste to cured cement.

Trial Kit #8716TK

(777-0560)Ea.

Contains: 1–3 Gm. syringe of translucent RelyX™ veneer cement, 1–2 Gm. syringe of translucent RelyX™ try-in paste, 20–0.12 mL unit doses of Scotchbond™ universal adhesive, 1–3 mL syringe of Scotchbond etchant, 25 etchant dispensing tips & 1 RelyX™ veneer shade guide.

Introductory Kit #8716

(777-0563)Ea.

Contains: 5–3 Gm. syringes of RelyX™ veneer cement (1 of each: A1, A3, white opaque, translucent & B0.5), 5–2 Gm. syringes of RelyX™ try-in pastes (1 of each: A1, A3, white opaque, translucent & B0.5), 1–5 mL vial of Scotchbond™ universal adhesive, 1–3 mL syringe of Scotchbond™ etchant, 48 disposable mix wells, 50 black microbrushes, 25 etchant dispensing tips & 1 RelyX™ veneer shade guide.

Refill Syringe #7614 3 Gm.

Specify:

Translucent (777-4162)

White Opaque (777-7126)

A1/Light Yellow (777-6771)

A3/Opaque/Yellow (777-7862)

B0.5/White (777-3623)

Try-In Paste

Syringe #7614T..... 2 Gm.

Specify:

Translucent (777-6699)

White Opaque (777-3740)

A1/Light Yellow (777-5570)

A3/Opaque/Yellow (777-2162)

B0.5/White (777-7833)

Rely X Ceramic Primer-#2721

(777-2015) 5 mL



SoloCem®

Coltene

Self-adhesive resin cement.

Formulated with MDP and 4-MET(A) monomers to provide high shear-bond and flexural strength without additional adhesive. Low shrinkage and water sorption. Indicated for ceramic, metal, and composite crowns/bridges and inlays/onlays, as well as endodontic posts and implant abutments. Contains zinc oxide. High radiopacity and color stability.

Introductory Kit

(380-0188)Ea.

Contains: 2–5 mL syringes (1 of each shade: dentin & translucent) & 20 brown mixing tips (10 of each size: short super fine & short fine).

Refill PackageEa.

Specify:

Dentin (380-0191)

White Opaque (380-0192)

Translucent (380-0194)

Contains: 1–5 mL syringe & 10 brown mixing tips (5 of each size: short super fine & short fine).



SpeedCEM® Plus

Ivoclar Vivadent

Self-adhesive, self-curing resin cement with light-curing option. Ideal for zirconia and metal ceramics. Bonds on moist and dry dentin in both self-cure and light-cure modes. Easy cleanup. High radiopacity.

System Kit

(945-0980)Ea.

Contains: 1-9 Gm. syringe (transparent), 1-5 Gm. Ivoclean bottle & 15 mixing tips.

Refill PackageEa.

Specify:

Transparent (945-0976)

Opaque (945-0977)

Yellow (945-0978)

Contains: 1-9 Gm. syringe & 15 mixing tips.



Variolink® Esthetic

Ivoclar Vivadent

Aesthetic light- and dual-curing luting composite for permanent cementation of ceramic and composite resin restorations. Patented, reactive, light-initiator Ivocerin® provides excellent shade stability, exceptional handling, and easy cleanup. Has lifelike fluorescence and is radiopaque.

VivaPen System Kits

Dual Cure

(945-0600)Ea.

Contains: 3-5 Gm. automix syringes (1 each: light, neutral & warm), 3-1.7 Gm. syringes of Try-In paste (1 of each: light, neutral & warm), 1-2 Gm. syringes of Total Etch, 1-2 mL Adhese Universal VivaPen, 1-5 Gm. bottle of Monobond Plus, 1-2.5 Gm. syringes of Liquid Strip, 30 VivaPen brush cannulas, 20 VivaPen protective sleeves, 20 Variolink® applicator tips, 1 brush holder & 50 soft brushes.

Light Cure

(945-0606)Ea.

Contains: 3-2 Gm. syringes (1 each: light, neutral & warm), 3-1.7 Gm. syringes of Try-In paste (1 of each: light, neutral & warm), 1-2 Gm. syringe of Total Etch, 1-2 mL Adhese Universal VivaPen, 1-5 Gm. bottle of Monobond Plus, 1-2.5 Gm. syringe of Liquid Strip, 30 VivaPen brush cannulas, 20 VivaPen protective sleeves, 20 Variolink® applicator tips, 1 brush holder & 50 soft brushes.

IPS e.max

Dual Cure

(945-0652)Ea.

Contains: 3-5 Gm. automix syringes (1 of each: light, neutral & warm), 3-1.7 Gm. syringes of Try-In paste (1 of each: light, neutral & warm), 1-2 Gm. syringe of Total Etch, 1-2 mL Adhese Universal VivaPen, 1-5 Gm. bottle of Monobond Etch & Prime, 1-2.5 Gm. syringe of Liquid Strip, 30 VivaPen brush cannulas, 20 VivaPen protective sleeves, 20 Variolink® applicator tips, 1 brush holder & 50 soft brushes.

Dual-Cure Refills

9 Gm.Ea.

Specify:

Light (945-0603) Warm (945-0605)

Neutral (945-0604)

Contains: 1-9 Gm. automix syringe & 15 automix tips.

Dual Cure Refills

5 Gm.Ea.

Specify:

Light + (945-0617) Warm (945-0620)

Light (945-0618) Warm + (945-0621)

Neutral (945-0619)

Contains: 1-5 Gm. automix syringe & 10 automix tips.

Light Cure RefillsEa.

Specify:

Light + (945-0611) Warm (945-0615)

Light (945-0613) Warm + (945-0616)

Neutral (945-0614)

Contains: 1-2 Gm. syringe & 5 applicator tips.

Try-In Paste RefillsEa.

Specify:

Light + (945-0622) Warm (945-0624)

Light (945-0623) Warm + (945-0625)

Neutral (945-0608)

Contains: 1-1.7 Gm. syringe & 5 applicator tips.

Invest in Your Practice...
NOW!
with the **Route 66 Lease Program**

- No payments for 6 months
- \$99 for the next 6 months
- Then level payments for 48 months*

Henry Schein Financial Services offers:

- Simple application process
- Flexible programs
- Affordable payments
- Fast approval process

To learn more, contact
Henry Schein Financial Services



1-800-268-6497 • hsf@henryschein.ca

*Minimum terms and conditions. Minimum transaction size: \$5,000. All programs feature a \$10 buy out. All transactions are subject to the satisfaction of underwriting guidelines, credit approval by third party lenders, and documentation requirements, and not all applicants will qualify. Certain other restrictions and additional terms and conditions apply. Offer expires December 31, 2021.

©2021 Henry Schein. Neither Henry Schein, Inc. nor Henry Schein Financial Services is a bank and neither conducts banking activities. Henry Schein Financial Services may receive a marketing fee from the vendor for products/services purchased.

STAY CONNECTED with



Stay informed on the most up-to-date news in the dental industry with our many social network platforms. Just one more reason to **Rely on Henry Schein.**

- www.Facebook.com/HenryScheinCanada
- [@hscheincanada](https://www.instagram.com/hscheincanada)
- [Henry Schein Canada](https://www.youtube.com/HenryScheinCanada)
- www.linkedin.com/company/HenryScheinCanada



C MAX 5 Ceramic Etching Gel

Gresco Products

Easy-to-use, yet effective etching gel. Used to etch both leucite and lithium-disilicate ceramics. Enhances bond strength by dissolving certain elements in the ceramic to create undercuts and voids in which the resin can flow.

(338-1802) 15cc



IPS Ceramic Etching Gel

Ivoclar Vivadent

Contains 4% hydrofluoric acid and designed to etch the internal aspects of all-ceramic IPS Empress CAD and IPS e.max CAD restorations. Prepares restoration for cementing.

Kit

(623-0160)Ea.

Contains: 1-5 mL bottle of IPS Ceramic Etching Gel & 1-30 Gm. bottle of Neutralizing Powder.

Refill

(623-5981) 5 mL



VITA® Ceramic Etch

VITA North America

This gel contains 5% hydrofluoric acid for etching ceramic restorations. Not for etching tooth structures.

Bottle

(809-0331) 6 mL



ONE COAT 7 UNIVERSAL

Coltene

Single-component, light-cured, universal adhesive. Bonds to both enamel and dentin. Application techniques include self-etch, selective etch, and total etch. Marginal seal helps minimize postoperative sensitivity. Works without silane or primer. Special One Coat 7.0 Activator ensures compatibility with chemical curing material.

Bottle Refill-5 mL

(889-0181)Ea.



Prime&Bond Elect®

Dentsply Sirona Restorative



Self-priming universal dental adhesive. Etchant, primer, and adhesive in a single bottle. Allows selective phosphoric acid etching of enamel only, phosphoric acid etching of both enamel and dentin, or no phosphoric acid etching prior to application. Outstanding bond strengths and protection against microleakage. Bonds resin-based materials to enamel and dentin as well as to metals and ceramic.

Introductory Kit

(222-0547)Ea.

Contains: 1-5 mL bottle of Prime&Bond Elect™, 1-3 mL syringe of tooth conditioner gel, 25 applicator tips, 50 flocked applicator tips & 1 dispensing well.

Dual Cure Introductory Kit

(222-0548)Ea.

Contains: 1-5 mL bottle of Prime&Bond Elect™, 1-3 mL syringe of tooth conditioner gel, 1-4.5 mL bottle of self-cure activator, 25 applicator tips, 50 flocked applicator tips & 1 dispensing well.

Bottle Refill

(222-0549)5 mL

Unit Dose 50-Pack

(222-0554)Ea.

Contains: 50 unit doses, 2 unit-dose holders & 50 flocked applicator tips.

Unit Dose 100-Pack

(222-0556)Ea.

Contains: 100 unit doses, 3 unit-dose holders & 100 flocked applicator tips.

OptiBond™ eXTRa Universal

Kerr Restoratives

Two-component, universal light-cure adhesive system. Self-etch, selective-etch, and total-etch techniques. Enriched with GPDM monomers, effectively etching the enamel while promoting the priming and sealing of dentin. Adhere to all composite, self-cure, dual-cure, and light-cure resin cement and core build-up materials. Indicated for direct and indirect restorations.

Bottle Introductory Kit

(133-6643)Ea.

Contains: 1-5 mL bottle of primer, 1-5 mL bottle of adhesive, 100 applicators & 25-4-cavity mixing wells.

Bottle Refill 5 mL

Specify:

Primer (133-6645) Adhesive (133-6642)

Unidose Introductory Kit

(133-6644)Ea.

Contains: 100-0.18 mL unidoses (50 of each: primer & adhesive) & 100 applicators.



Optibond™ Universal

Kerr Restoratives

Single-component, light-cure adhesive. Provides adhesion to a variety of surfaces and substrates for both direct and indirect applications. Can be used in self-etch, selective-etch, and total-etch techniques.

Bottle Kit

(127-5562)Ea.

Contains: 1-5 mL bottle of Optibond Universal, 1-3 Gm. gel-etchant syringe, 10 gel-etchant syringe tips, 25 mixing wells & 100 applicator tips.

Bottle Refill

(127-5564) 5 mL

Unidose® Kit

(127-5565)Ea.

Contains: 100-0.18 mL unidoses, (18 mL total) 1-3 Gm. gel-etchant syringe, 20 gel-etchant syringe tips, 25 mixing wells & 100 applicator tips.

Gel Etchant Complete Kit

(123-8669)Ea.

Contains: 3-3 Gm. syringes & 30 disposable dispensing tips.



3M

Scotchbond™ Universal Adhesive

A combined total-etch, self-etch, and selective-etch adhesive. Consistent bond strength to both moist and dry etched dentin. Single-bottle system for both direct and indirect restorations. Provides consistent, reliable bond strengths with virtually no postoperative sensitivity.

Vial Delivery Systems

Introductory Kit #41254

(378-0348)Ea.

Contains: 1-5 mL vial of universal adhesive, 1-3 mL syringe of etchant, 50 mixing wells, 50 microbrushes & 25 dispensing tips.

Refill Pack-#41258

(378-0349) 5 mL

Unit-Dose Delivery Systems

Refill Pack #41255

(378-0352)Ea.

Contains: 40-0.12 mL unit doses, 1-3 mL syringe etchant & 25 dispensing tips.

Giant Pack #41256

(378-0354)Ea.

Contains: 100-0.12 mL unit doses, 1-3 mL syringe etchant & 25 dispensing tips.

Bulk Pack #41257

(378-0353)Ea.

Contains: 200-0.12 mL unit doses, 1-3 mL syringe etchant & 25 dispensing tips.

Dual-Cure Activator-#41264

(378-0350) 5 mL

Etching Gel Package #41263

(777-0453)Ea.

Contains: 2-3 mL syringes etchant & 50 dispensing tips.





SMART MATERIAL FOR SMART APPLICATIONS

**EXCLUSIVE
OFFER**
Buy 1, **Get 1**
**Optiglaze
Color
(Clear)***

CERASMART® 270

Force Absorbing Hybrid Ceramic CAD/CAM Block

CERASMART® 270 is an aesthetic hybrid block made to provide fast and long-lasting indirect restorations.

- Mills quickly and gentle on milling burs
- Virtually no chipping
- No need for sintering/crystallization firing
- Flexible, tough and fracture resistant
- Easy to polish, retains gloss over time and is self polishing
- Highly wear resistant

More Shades Available!

CERASMART® 270 for CEREC® 5-Block Package (Size 12): \$149.99 4@ 145.99

(333-5258) HT/A1 (333-5260) HT/A2 (333-5261) HT/A3 (333-5220) LT/A1

CERASMART® 270 for CEREC® 5-Block Package (Size 14): \$159.99 4@ 155.99

(333-5272) HT/A1 (333-5275) HT/A2 (333-5277) HT/A3 (333-5281) LT/A1

CERASMART® 270 for CEREC® 5-Block Package (Size 14L): \$177.99 4@ 172.99

(333-5246) HT/A1 (333-5247) HT/A2 (333-5248) HT/A3 (333-5251) LT/A1



(333-5267) LT/A2

(333-5268) LT/A3

(333-5284) LT/A2

(333-5285) LT/A3

(333-5252) LT/A2

(333-5254) LT/A3

www.gcamerica.com • www.gcamerica.com/training • © 2020 GC America Inc.

*GC America code 2471-2. To receive your no charge goods, send a copy of your Henry Schein Dental invoice dated 12/27/20 - 3/27/21 to GC America, 3737 West 127th Street, Alsip, IL 60803 or fax to 800-423-2963 Attn: Customer Care Group. Must be redeemed within 6 weeks of Henry Schein Dental invoice date. Offers cannot be combined with any other offers or promotions. No substitutions allowed. Offers are valid in the US and Canada only. No charge goods will be of equal or lesser value to the lowest selling unit price. Please allow 6-8 weeks for no charge goods delivery. No charge goods will be fulfilled by GC America unless noted otherwise. All promos expire 3/27/21. †CEREC® is a registered trademark of Dentsply Sirona. PlanMill® is a registered trademark of Planmeca Oy.



We Know Digital Production



keyprint®

KeySplint Soft®



keymill™

Denture Base Disc



keymill™

Wax Disc



enamelite®
by keystone industries

CAD/CAM High Glaze and Scanning Spray



Finishing



Polishing



Henry Schein Cares

Helping Health Happen®

MAKING A DIFFERENCE

Together, we are Helping Health Happen®



Back to School

Providing hundreds of backpacks filled with school supplies and clothing annually to give underprivileged children a sense of confidence as they head back to school.



Community

Helping underserved people in communities across Canada receive dental care they normally could not afford.



Disaster Relief

Shipping more than \$1 million of healthcare products a year to aid when natural disasters strike globally.



Going Green

Reducing paper, planting trees and saving energy. Our Eastern Distribution Centre has gone paperless, eliminating 17,000 lbs of paper usage each year!



Henry Schein Cares

Helping Health Happen®



1-800-668-5558 | www.henryschein.ca



IMPRESSION MATERIAL



Blu-Bite® for CAD/CAM



Vinyl polysiloxane material specifically designed to optimize scanning and image processing using CAD/CAM restoration systems. White material produces greater reflectibility than colored material. Achieves a rigid set while maintaining resistance to fracture. No contrast spray required. Superfast setting (<30 seconds).

Complete Package Ea.
(112-5483)

Contains: 2-50 mL cartridges & 6 pink mixing tips.



Virtual® CADbite

Ivoclar Vivadent

Scannable, vinyl-polysiloxane impression material. Used for the optical recording of data when fabricating restorations with a CAD/CAM system. Nonslumping consistency, short setting time, and good dimensional stability. High final hardness of 32 Shore-D and high fracture toughness. Mint scent.

Refill Pack Ea.
(945-9137)

Contains: 2-50 mL cartridges & 12 mixing tips.



Flexitime® Fast & Scan

Kulzer

Scannable vinyl polysiloxane material that supports CAD/CAM manufacturing of prosthetic restorations and offers optimal workflow integration for 1- to 3-unit impressions and their digitalization. Increased efficiency in the handling process with the laboratory.

50 mL Cartridge Refill Ea.

Specify:

Light Flow (171-0496)

Medium Flow (171-0501)

Contains: 2-50 mL cartridges & 12 mixing tips.

Easy Putty Refill Package

(171-0487) Ea.

Contains: 2-300 mL cartridges (base & catalyst) & 1 scoop.

380 mL Dynamix® Refill Ea.

Specify:

380 mL Dynamix® Refill

Heavy Tray (171-0489)

Contains: 2-380 mL cartridges, 2 Dynamix mixing tips & 1 fixation ring.

CAD/CAM ACCESSORIES

MILLING TOOLS



CEREC® Milling Instruments

Dentsply Sirona

For Use with CEREC® and InLab® Milling burs that do the actual milling of the composite blocks used in the milling machine. For Sirona MC, MC X, Practice Lab, and inLab® MC XL.

12mm 6 per Pkg. **205.99**
(811-1052)

HENRY SCHEIN®
TECHNICAL SERVICES
1-888-818-5838



New!

CEREC® Primemill Milling Instruments

Dentsply Sirona

Designed for grinding of glass ceramic and composite materials. Must be used with the appropriate torque wrench.

Diamond CS..... 6 per Pkg.

Specify:

1.4mm (137-6130) 1.0mm (137-6132)

1.2mm (137-6131) 0.6mm (137-6137)

CEREC® Primemill Burs

Dentsply Sirona

Designed for milling of zirconia and PMMA materials. Must be used with the appropriate torque wrench. Can be used for wet or dry milling.

ZrO₂ CS-2.5mm 3 per Pkg.
(137-6138)

CS-1.0mm 3 per Pkg.
(137-6139)

CS-0.5mm 3 per Pkg.
(137-6140)

PrograMill® Milling Tools

Ivoclar Vivadent

For use with PrograMill® One Specifically designed for different CAD/CAM materials. Innovative color concept with color-coded rings on the shaft. Color coding is also used in the PrograMill One App. Allows intuitive and easy tool management.



Milling Bur for Glass Ceramics

Red Ring..... Ea.

Specify:

1.2mm (945-2206) 0.5mm (945-2207)

Rely on Us™

Sirona MC X5 Milling Burs

Dentsply Sirona

For Use with inLab® MCX5

Burs used for the milling and grinding process for workpieces in block or round form. Grinding burs are used in the wet milling process; milling burs can be used in the wet or dry milling process.



Milling Bur for Zirconia Oxide (ZrO₂) Diamond Coated

Yellow BandEa.

Specify:

2.5mm (811-6943) 0.5mm (811-6937)

1mm (811-6940)



Milling Bur for Composite Coated-Blue Band

2.5mmEa.

1mmEa.

0.5mmEa.



Milling Bur for Sintering Metal-Gray Band

2.5mmEa.

1mmEa.

0.5mmEa.



Milling Burs for PMMA Plastic-Red Band

2.5mmEa.

1mmEa.

0.5mmEa.



Grinding Tools for Glass Ceramics

White BandEa.

Specify:

2.2mm (811-6911) 1.2mm (811-6915)

1.4mm (811-6913)

0.6mmEa.



Two Striper® Milling Tools

Premier Dental

Milling tools for use with E4D® Mill/PlanMill 40. Feature the same faster, cooler cutting, and longer lasting qualities that have made Two Striper diamond instruments the market leader. They are manufactured using the patented P.B.S.® bonding process that maximizes the exposure of the diamond-cutting surfaces. Natural, virgin diamond crystals on Two Striper rotary instruments contain more corners and angles than synthetic diamonds used on electroplated tools.

Burs-Standard Length ...6 per Box

Specify:

Tapered (378-6546)

Ellipsoidal (378-0560)

Conical (378-1031)

Burs-Long Length6 per Box

Specify:

Tapered (378-0867)

Conical (378-0868)

Assorted Kit-Standard Length

(378-0206)Ea.

Contains: 6 milling tools (2 of each shape: tapered, ellipsoidal & conical).

SCANNING & PREPARING

CAD/CAM Scanning Spray



Covers material with a thin, uniform glaze to allow for a more accurate scan and therefore, a more accurate fit. Spray prior to scanning and only use for nonintraoral scanning, such as on crowns, models, and impressions. Extension spray nozzle creates a tighter spray pattern.

Complete Package

(570-0541)Ea.

Contains: 1 of each: 3.1 oz. spray can & nozzle.

CEREC® DENTATEC Cleaner & Lubricant

Dentsply Sirona

Concentrated milling solution used to clean and lubricate CEREC® and inLab® systems. Dilute with a ratio of water, depending on the unit and desired effect.

Bottle

(811-0775) 1000 mL



CEREC® Optispray

Dentsply Sirona

Enables easy and perfect optical impressions of CAD/CAM restorations and aids in the topographical recordings of prepared and unprepared teeth and their surroundings. Thin, homogeneous contrast layer helps deliver precise digital-impression results: simply coat the tooth and trigger the camera. Spray head can be turned or angled to the desired position, depending on application. Blue color clearly indicates sprayed area.

50 mL Spray

(587-2706)Ea.

Contains: 1 of each: 50 mL can, special nozzle & stabilization tube.

200 mL Spray

(587-2751)Ea.

Contains: 1-200 mL can, 3 special nozzles & 1 stabilization tube.



CEREC® Stone BC

Dentsply Sirona

Super-hard stone material that allows for highly precise, scannable models. Features optical properties such as brightness and contrast. Use in combination with the CEREC® Bluecam.

1.2 Kg.

(587-2713)2 per Pkg.

0.1 Kg.

(587-2712)20 per Pkg.



Composite Block Defoaming Solution

D4D Technologies

Specially formulated for PlanMill 40 when using composite blocks, to alleviate possible foaming. Applied into the chamber as needed.

4 oz. Bottle

(681-4146)Ea.



Digital Impression Accessories

3.75 Gm. Powder for Chairside Oral Scanner (777-6082)Ea.

Powder Sprayer Tips for Chairside Oral Scanner (777-0253)3 per Pkg.

Powder Sprayer for Chairside Oral Scanner (777-0727)Ea.

Contains: 1 sprayer, 6 tips, 1 battery, 1 screwdriver & 2 brushes.





Earth Stone™ Dental Stone

Tak Systems

Gypsum composite that takes only 1 min. to prepare and 10 min. to bring to full working hardness. "Shake & Pour" feature eliminates dust hazards. No bowls to wash or mess to clean. Bubble-free models can be prepared without a vibrator at chairside.

Complete Package

(966-2932) 18 per Box
Contains: 18 of each: 90 Gm. envelopes & 30 mL activator liquids.



Enamelite™ CAD/CAM Scanning Spray

Keystone Industries

CAD/CAM scanning spray that covers materials with the thinnest, most uniform and accurate layer. This allows for a more accurate scan, and in turn, a more accurate fit. It is meant to be sprayed prior to scanning, and should only be used for nonintraoral scanning, such as on crowns, models, and impressions.

Spray Can

(242-0946) 3.1 oz.

Planmeca Emerald™ Accessories

Planmeca

Tip

(631-0182) Ea.

Slim Line Tip

(631-0287) Ea.



PlanScan Scanner Tips

Planmeca

Removable tips for the Planscan system. Standard tips are ideal for general scanning procedures and a larger field of view. Small portrait tips are ideal for patients with smaller mouths, for enhanced reach to the posterior, for compromised dentition, or for capturing interproximal areas. Small landscape tips are ideal for anterior scanning and easy capture of buccal bite. Tips require effective high-level disinfection.

Standard-Color

(631-0056) Ea.

Standard

(631-5426) Ea.

Autoclavable Ea.

Specify:

Standard Size (445-7436)

Small Landscape (445-4531)

Small Portrait (445-2695)

PrograMill® Fluid

Ivoclar Vivadent

Concentrated milling coolant.

Specifically formulated for the PrograMill One milling machine.

(945-2209) 1000 mL



Scan-Mate™ Retraction Instrument

Armor Dental

Nonreflective retraction instrument for intraoral scanning and any retraction need. The soft silicone tip provides more comfort to the patient, and the flexible neck can be bent to provide complete access to all areas of the mouth, reducing scanning time for crowns, bridges, implants, orthodontics, and many other procedures. Autoclavable.

(131-7726) 4 per Box

Scanspray

Renfert USA

High-end surface matting that provides a very thin layer application and has uniform grain size. Special valve technology reduces clogging.

(228-2324) 200 mL



TRIOS® Scanner Tips

3Shape

For use with the TRIOS® Standard and Color Digital Impression Systems.

Autoclavable 10 per Box

Specify:

TRIOS® Color (820-2723)

TRIOS® Standard (820-4619)


Scanner Tips with Mirror

(820-0070) 3 per Box

Scanner Protection Tip

(820-4055) Ea.

MISCELLANEOUS



Glass Mixing Slab

HENRY SCHEIN®

#4, 4" x 2½" x ¼" Thick

(100-4432) Ea.

6" x 3" x ½" - Thick #H

Polished on both sides. Beveled edges.

(100-4129) Ea.

6" x 3" x ¾" - Thick #T

Polished on both sides. Beveled edges.

(100-3163) Ea.



Primescan Scanner Accessories

Dentsply Sirona

Autoclavable Sleeve

(939-7249) Ea.

Disposable Window

(939-7250) 50/Pk.



CONTACT YOUR HENRY SCHEIN® REPRESENTATIVE FOR ALL YOUR

CAD/CAM MATERIAL NEEDS

HENRY SCHEIN®
SOLUTIONS FOR HEALTH CARE PROFESSIONALS

Rely on Us

1-800-668-5558 • www.henryschein.ca



Sable Brushes



- Used for porcelain and staining
- Lightweight: reduces hand fatigue

StraightEa.

Specify:

000 (100-7603)	1	(100-1190)
00 (100-4238)	2	(100-3793)
0 (100-3510)		

Size 3
(101-8789) Ea.

Size 4
(101-8821) Ea.

Size 6
(101-8859) Ea.

Angled Ea.

Specify:

000 (100-2092)	0	(100-7731)
00 (100-6668)		

Size 1
(100-8149) Ea.

Note: Sizes are in ascending order. 000 is the smallest.



Ivoclar Vivadent Block Locker

Ivoclar Vivadent Storage unit for CAD/CAM blocks and accessories. Allows for customized organization and easy access to materials. For use with both Dentsply Sirona CEREC® and Planmeca PlanMill™ materials.

Small Cabinet-19.5" x 16" x 9.25"
6 internal organizers and 6 block dividers. Holds up to 288 individual blocks.

(945-2067) Ea.

Large Cabinet-25" x 26.6" x 9.25"
12 internal organizers and 4 block dividers. Holds up to 48 packages of blocks or 312 individual blocks.

(945-2065) Ea.



Pic-n-Stic™

Pulpdent
2½" long plastic stick with an adhesive tip at one end. Acts as a handle for small objects. Can be bent 35°-45°. Great for picking up inlays, veneers, bands, brackets, files, reamers, and other items.

(384-3960) 60 per Box
Contains: 4 packs of 15 sticks.



CAD/CAM Block Locker

Zirc

Storage unit for milling blocks and accessories. Wall-mounted unit with a gray, powder-coated finish. Features internal organizers with 7 different color options to fit interior design of practice. Drawer accommodates additional accessories. Recommended for most block sizes.

Small Cabinet-19½" x 16" x 9¼"
6 organizer trays hold up to 288 blocks.
Small..... Ea.

Specify:

Copper	(658-0451)
Midnight Blue	(658-0448)
Vibrant Blue	(658-0453)
Vibrant Orange	(658-0454)
Vibrant Pink	(658-0452)
Plum	(658-0450)
Teal	(658-0449)

Large Cabinet-25" x 26⅝" x 9¼"
12 organizer trays hold up to 48 packages of blocks or up to 312 blocks.

Large..... Ea.

Specify:

Copper	(658-0458)
Midnight Blue	(658-0455)
Vibrant Blue	(658-0460)
Vibrant Orange	(658-0461)
Vibrant Pink	(658-0459)
Plum	(658-0457)
Teal	(658-0456)



CAD/CAM Block Organizer™

Zirc

Organizer trays to store packages of blocks or individual blocks utilizing existing drawer space. Organizer..... Ea.

Specify:

Copper	(658-0484)
Midnight Blue	(658-0483)
Vibrant Blue	(658-0480)
Vibrant Orange	(658-0481)
Vibrant Pink	(658-0482)
Plum	(658-0478)
Teal	(658-0479)



Flex Tab™ Flexible Clearance Tabs

Kerr Lab

Easy and accurate way to make direct intraoral measurements of occlusal and interproximal tooth-preparation clearances. Place an individual tab between a prepared tooth and its opposing occlusal surface to ensure the laboratory has sufficient clearance to develop the perfect restoration. Tabs are color coded by thickness.

Individual Sizes..... 30 per Roll

Specify:

Pink-1mm	(236-4843)
Green-1.5mm	(236-4871)
Blue-2mm	(236-4923)

Assorted Package

(236-4324) 90 per Pkg.
Contains: 30 tabs of each: 1mm, 1.5mm & 2mm.

Clearance Guide

(236-1660) 2 per Pkg.



OptraStick®

Ivoclar Vivadent

Flexible plastic placement instrument. Features a ball-shaped, flexible adhesive tip. Tip allows restorations and small objects to be picked up, held, carried, and placed quickly and easily, while the thin part of the stick can be easily bent and flexed in all directions.

(945-3227) 50 per Pkg.



New!

Replacement Filter

Ivoclar Vivadent

Replacement filter for PrograMill® One milling machine.

(945-2208) Ea.



Practice Care is Our Priority

At Henry Schein Canada, our mission is to improve the lives of those we touch by focusing on practice care, so dental professionals can focus on patient care.

Practice Care is a combination of efficient solutions and integrated technologies designed to help you operate a productive practice, attain your business goals, and assist in the delivery of quality patient care.

Our Solutions focus on connecting all of your practice needs for a streamlined, digital workflow.



The way we put it all together helps you operate a productive practice and deliver the best possible patient experience.

To place your order, contact your
Henry Schein Sales Consultant or call:

1-800-668-5558

Order online 24/7 at:

henryschein.ca



Henry Schein Canada




Save Money. Save Time.

Mill Custom Titanium
Abutments In-house.



AK-1 Abutment Kit For the DWX-42W Wet Dental Mill

Milling custom titanium abutments is now easier and more cost-effective than ever, with a DWX-42W dental mill and the AK-1 Titanium Abutment Kit, from DGS SHAPE.



1
Replace regular DWX-42W tank tray with AK-1 kit tray.





2
Connect AK-1 material adapter to material holder in the DWX-42W.



3
Insert and attach a titanium block, then begin milling.

See how easy and affordable it is to add custom titanium abutment milling to your lab output at www.rolanddga.com/AK-1

Get the Planmeca FIT®
Complete CAD/CAM Solution
as low as \$112,000* CAD



SCAN CODE
TO LEARN MORE



Affordable
same-day dentistry
is here!

PLANMECA

© 2021 Planmeca U.S.A. Inc. All Rights Reserved. PMCA-CDM-AD0003A 03/21
*Prices are listed in Canadian dollars and subject to change based on current exchange rates.

Rely on Us™

HAVE YOU
SWITCHED TO
**HENRY SCHEIN
BRAND**
PRODUCTS YET?



Over 5,000 **Henry Schein Brand** products are the perfect choice for exceptional quality and maximum value.

**Now, more than ever, you can rely on
HENRY SCHEIN BRAND products:**

- **Pay less** and **save more**
- Made by **leading manufacturers**
- **Quality** that meets or exceeds competing products
- Backed by our **100% satisfaction guarantee** – try them today with no risk!



Why pay more when you don't have to?
Contact us today to find out how much you can save!

1-800-668-5558 | henryschein.ca