ULTICLEAN 500L

Dental Scaler + Air Polisher



Modern approach for prevention, cleaning and minimally invasive prophylaxis treatments.



SIMPLE & SOPHISTICATED

Ease of use, a complete unit that comes with all the tools to allow for many procedures:

AIR POLISHING + PERIODONTAL TREATMENT + ENDODONTIC TREATMENT + IMPLANT MAINTENANCE + EFFICIENT AND COMFORTABLE PLAQUE REMOVAL

AUTOMATIC IDENTIFICATION OF PROCEDURE MODE

Automatically recognizes which feature you require when you pick up either handpiece – scaling or air polishing.

AIR POLISHING

Compressed air precisely delivers powder and water effectively removing surface substances and polishes the tooth surface. It is less harsh on the tooth enamel with minimal abrasion. Delivers an efficient and comfortable treatment ensuring a pleasant experience for patients for supragingival and subgingival treatments.

One package with everything you need to have clinical confidence.
Two-year warranty.



(940-2173) UC500L

(940-2580) Cart for UC500L (sold separately)

one

FROM PERFORMANCE TO PEACE OF MIND... package

B.A.INTERNATIONAL

Certified to EN ISO 13485:2016 Quality Management System for Medical Devices

Exclusively available from

HENRY SCHEIN®

www.henryschein.ca | 1-800-668-5558

B.A. INTERNATIONAL
Certified to EN ISO 13485:2016, Quality Management System for Medical Devices



HENRY SCHEIN®

Hygiene: Tips, Accessories, Powders



Titanium alloy tip

The hardness of titanium alloy tip is lower than that of cementum and enamel, so teeth are not damaged.

VIEW	ITEM CODE	HENRY SCHEIN ORDER NUMBER	SIZE	PROCEDURE TYPE	DESCRIPTION
	BA150180	(940-2735)	BAC33E	Periodontal	Remove calculus deep in the periodontal pocket and rinse the periodontal pocket.
	BA150181	(940-2737)	BAC50RE	Periodontal	For subgingival calculus detection and subgingival scaling.
	BA150182	(940-3182)	BAC50LE	Periodontal	For subgingival calculus detection and subgingival scaling
	BA150184	(940-3180)	BAC59E	Periodontal	For subgingival calculus detection and subgingival scaling.
	BA150183	(940-3178)	BAC52E	Scaling	For supragingival scaling and removal of dental calculus / plaque.
	BA150185	(940-3179)	BAC56E	Scaling	For supragingival removal of larger calculus.
	BA150195	(940-2747)	BAC60E	Endodontic	For preparation and irrigation of root canal.
	BA150196	(940-2745)	BAC62E	Endodontic	For Root canal irrigation.
	BA150188	(940-2736)	BAC90E	Implant	For cleaning implant and restoration.
	BA150186	(940-2744)	BAC94E	Implant	For cleaning implant and restoration.
	BA150187	(940-2742)	BAC95E	Implant	For cleaning implant and restoration.

Accessories



ITEM CODE	
(940-2758)	BAC5P Air Polisher Handpiece for supragingival treatment
(940-2759)	BAC5B Air Polisher Handpiece for subgingival treatment
(940-2760)	BAC5L Ultrasonic Handpiece scaler
(940-2757)	Single use Nozzles (pack of 10) (Subgingival Air Polisher Tips)
(940-2580)	Cart



Ulticlean Air Polishing Powders

- NEW polishing powders, available in packs of 4 in lemon, mint and perio plus
- Mint & Lemon powders: 40µm, sodium bicarbonate base



(940-0636) Unflavoured
(940-0638) Mint
(940-0637) Lemon

Exclusively available from HENRY SCHEIN® www.henryschein.ca | 1-800-668-5558