ULTICLEAN 500L Dental Scaler + Air Polisher

WHAT COMES IN THE BOX

INCLUDED IN THE PURCHASE OF THE UC500L

- 2 x LED Ultrasonic handpieces
- 1 x Air Polishing handpiece
- 1 x Perio Air polishing handpiece
- 2 x Water tanks (1500 ml & 600 ml)
- 2 x Powder tanks (1 Sub & 1 Supra)
- 14 x Tips
    - 2 x Subg scaling (59E)
    - 2 x Heavier Calc removal (56E)
    - 2 x Root canal (60E)
    - 2 x Sub G scaling right (50RE)
    - 2 x Sub G scaling left (50LE)
    - 2 x Implant restoration (94E)
    - 2 x Calc removal within perio pockets (33E)
- 1 Box of 10 single use perio tips
- 1 x Foot pedal
- 1 x Torque wrench
- 3x Sterilization boxes (2 x scalers, 1 x air polisher)
- 2 x Cleaning needles (1 small/1 large)
- 4 x Compressor connectors (2 female, 1 male, T-junction)
- 1 x Filter wrench
- 1 x Replacement LED module
- 1 x Replacement Oring kits

ULTICLEAN 500L

Dental Scaler + Air Polisher

Titanium alloy tip



VIEW	ITEM CODE	SIZE	PROCEDURE TYPE	DESCRIPTION
	BA150180	BAC33E	Periodontal	Remove calculus deep in the periodontal pocket and rinse the periodontal pocket.
	BA150181	BAC50RE	Periodontal	For subgingival calculus detection and subgingival scaling.
	BA150182	BAC50LE	Periodontal	For subgingival calculus detection and subgingival scaling
	BA150184	BAC59E	Periodontal	For subgingival calculus detection and subgingival scaling.
	BA150183	BAC52E	Scaling	For supragingival scaling and removal of dental calculus / plaque.
	BA150185	BAC56E	Scaling	For supragingival removal of larger calculus.
	BA150195	BAC60E	Endodontic	For preparation and irrigation of root canal.
	BA150196	BAC62E	Endodontic	For Root canal irrigation.
	BA150188	BAC90E	Implant	For cleaning implant and restoration.
	BA150186	BAC94E	Implant	For cleaning implant and restoration
	BA150187	BAC95E	Implant	For cleaning implant and restoration

- BAC59E is the traditional Piezo Tip used in 90% of all scaling procedures sub & supra
- BAC56E is larger in diameter than 59E and can help remove supragingival heavy calculus with a ball tip for greater impact
- BAC33E is a finer tip used for deep perio pockets