
Essentials

FOR ALL YOUR PRACTICE NEEDS

Compact LED Curing Light



- Ideal for the hygiene operatory
- Easy-to-use, cordless, lightweight design
- 1000 mW/cm^2 light intensity
- Full power curing times of 5, 10, & 15 seconds
- 8 mm black-coated light guide
- Universal micro-USB charger
- Audible end-of-cycle notification tone
- 1-year warranty

COMPACT LED CURING LIGHT
Contains: 1 Essentials Compact LED Curing Light, 1 Holding Cradle, 1 Micro-USB Cable, 1 Power Adapter, 1 Glare Shield, 2 Anti-glare Cups

(570-0571). . . . . Ea.

HENRY SCHEIN®

1.800.668.5558
www.henryschein.ca