Essentials

CORE SUPPLIES. EVERYDAY LOW PRICES.

SECURELY SEALED. STERILIZATION MAINTAINED.

NEW, with dual indicators



ALSO AVAILABLE IN NEW SIZES!

ESSENTIALS SELF SEAL STERILIZATION POUCHES

WITH DUAL INDICATORS

FIND THESE PRODUCTS AND MORE AT:

HENRYSCHEIN.CA

Essentials

CORE SUPPLIES. EVERYDAY LOW PRICES.

Essentials Self Seal Sterilization Pouches

WITH DUAL INDICATORS



- Pouches are mechanically sealed using proprietary technology that fuses green-tinted film with high-grade medical paper, creating a strong bond
- Dual indicators change color when a sterilization process has begun and to confirm sterilization use a biological monitor
  - Steam Autoclave: Turns from pink to cocoa at 121°C (air-removal by gravity) after a minimum of 30 minutes with 25 minutes drying time or 134°C (air-removal by Pre-Vacuum) after a minimum of 4 minutes with 20 minutes drying time.
  - Ethylene Oxide (EtO): Turns from blue to golden brown after 120 minutes at 50°C, RH 60-85%, and EtO concentration of 600 mg/L. Aeration at 25°C: Single pouch: 5 days. Double pouches 3 days.
- Reliably meets sterilization standards when used according to the Instructions for Use (IFU)
- Meets ISO 11607 and ISO 11140-1 standards
- Made with leak resistant, medical grade paper with easy to use thumb notch
- When properly stored, instruments maintain their sterility for up to 1 year post sterilization
- Compatible with autoclaves and ethylene oxide (EtO) chamber sterilizers

SIZE	QUANTITY	ITEM CODE
2¼" x 4"	200/Box	(570-4241)
3½" x 5¼"	200/Box	(570-4240)
2¾" x 9"	200/Box	(570-4238)
3½" x 9"	200/Box	(570-4236)
3½" x 9"	500/Box	(570-4244)

	SIZE	QUANTITY	ITEM CODE
	5¼" x 10"	200/Box	(570-4237)
	4¼" x 11"	200/Box	(570-4239)
NEW!	**7½" x 13"**	200/Box	(570-4243)
NEW!	**9" x 15"**	100/Box	(570-4242)

ORDER TODAY AT HENRYSCHEIN.CA